
AlaxoLito

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
Gebrauchsanweisung
Instructions for Use
Manuel d´utilisation
Istruzioni per l´uso
Gebruiksaanwijzing

AlaxoLito

Nasenstent

+PLUS

**Bestellnummer der AlaxoLito Produktvarianten:
Order numbers for AlaxoLito product variants:
Numéros de commande des variantes de l'AlaxoLito:
Numero d'ordine delle varianti del prodotto AlaxoLito:
Bestelnummers van de AlaxoLito-productvarianten:**

REF

AL-002-003 AlaxoLito Set
AL-002-006 AlaxoLito Plus 70 Set

Deutsch	3
English	10
Français.....	16
Italiano.....	23
Nederlands	30

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1 Instructions

The AlaxoLito Nasal Stent is a CE marked prescription medical device. It must only be used according to the briefing of the attending doctor.

1.1 CE Mark

The product bears the CE mark. The declaration of conformity assures that the product complies with the requirements of the EU directive for medical devices 93/42/EEG. According to the classification rules in Annex IX of the EU directive the AlaxoStent is a Class I product.

1.2 General Instructions

- The instructions for use is a component of the product. It has to be accessible to the user near to the product. Close adherence to these instructions is precondition for intended use and correct handling of the product.
- In case of damage resulting from use of third-party components or consumables no warranty claims apply.
- Alaxo is only liable for the product in respect of safety, reliability and functionality if it is used in compliance with the instructions for use.
- The instructions for use corresponds to the design of the product and the applicable safety norms as of going to print. All rights are reserved for the specified parts and components, processes, names and products.
- Full or partial reproduction subject to prior written approval by Alaxo.

1.3 General safety instructions

- The product has been developed and designed by Alaxo in such a manner that hazards during intended use are eliminated as far as possible. Nevertheless, we strongly recommend to strictly adhere to the following safety precautions in order to further reduce any residual risks to a minimum.
- It is not allowed to modify or reconstruct the product. In that case the CE mark lapses. Use of altered products is prohibited. In addition, Alaxo cannot take any warranty or liability for modified or reconstructed products. In order to assure safe application of the product users are liable for adherence to the instructions.
- The original packaging should be stored for any potential return shipment. Only the original packaging ensures optimal protection of the product during transport. Should return shipment become necessary during the warranty period Alaxo does not take any liability for damages resulting from insufficient packaging.
- The product must not be accessible to children.
- The user has to inspect the product for functional safety and good condition each time prior to application.
- The user has to be familiar with handling of the product.
- The cleaning instructions have to be observed and obeyed to. Other cleaning procedures must not be applied. The cleaning spray (Prontosan® Wound Spray) must only be used for disinfectant cleaning of the AlaxoLito Nasal Stent and the accessories. Do not swallow!

1.4 Symbols and indications

- One Prontosan[®] Wound Spray bottle may be used for 12 months after first use.
- The AlaxoLito Nasal Stent may be used for 18 months after first use.
- In the instructions for use, on packagings and on the product the following figures and symbols are used:

Consult instructions for use

Batch code

Month/Year Use by: Month/Year

Order number

Manufacturer

CE Mark

Temperature limitation: store between 4° and 50°C

Store inaccessible for children

1.5 Precautionary measures

The attending doctor has to explain handling of the AlaxoLito Nasal Stent to the patient. Especially, the doctor has to explain correct application and positioning of the AlaxoLito stent in the nose and its removal.

The attending doctor has to advise on potential risks and complications connected to the AlaxoLito Nasal Stent therapy and has to evaluate if the AlaxoLito therapy is applicable to the specific patient.

2 Product information

2.1 Intended use

- The AlaxoLito Nasal Stent is a medical device and self-applied by the patient. The AlaxoLito stent is positioned in the anterior nasal cavity in order to mechanically splint the anterior nasal airway.
- The AlaxoLito Nasal Stent serves to prevent nasal valve and/or nasal alar collapse.
- The AlaxoLito Nasal Stent splints the anterior inferior or middle nasal passage. The AlaxoLito Plus Nasal Stent splints the inferior or middle nasal passage over its entire length.
- The AlaxoLito Nasal Stent can be used to treat nasal valve or nasal alar collapse and turbinate hyperplasia turbinate hyperplasia and potentially resulting snoring.
- The AlaxoLito Nasal Stent is a CE marked and prescription medical device.
- The attending doctor has to examine if the AlaxoLito Nasal Stent is an appropriate therapy for the patient and if the nasal cavities and passages of the patient are appropriate for application of the AlaxoLito stent.
- The AlaxoLito Nasal Stent must be used only by one single patient.

- The AlaxoLito Nasal Stent is suitable for repeated use up to 18 hours per application.
- The AlaxoLito Nasal Stent must be used only by patients of age of 18 years and higher.
- The AlaxoLito Nasal Stent is not suitable for life-supporting application in emergency care.

2.2 Contraindications

- Dry rhinitis
- Strong deviations of the nasal septum
- Disposition for nose hemorrhage
- The AlaxoLito Nasal Stent shall not be used by persons with reduced mental ability, with reduced motor skills, with psychic restrictions, mental disorders or spastic syndromes.

2.3 Warnings

- Persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this stent.
- Persons with allergic reactions to polyhexanide or undecylenamidopropylbetaine may suffer an allergic response to the cleaning spray.
- Do not use the AlaxoLito Nasal Stent if this leads to pain! Do not terminate therapy without consulting your attending doctor.

2.4 Unintended use

- The only permitted use is splinting of the airway in the nasal cavity. Other uses are prohibited.
- The AlaxoLito Nasal Stent must not be introduced into orifices of the body other than the nose (e.g. ears or mouth).
- The AlaxoLito stent must not be pulled over any body parts or limbs.
- The product may only be used for the approved medical purposes. It must not be used for any other purposes, especially non-medical purposes.

2.5 Product description

The AlaxoLito Set comprises the following components (see **Fig. 1**):

- 1) Self-expanding nitinol braid [1] with clip lock [1a] (nitinol is a shape-memory metal, a nickel-titanium alloy widely used for medical devices)
- 2) Pusher [2] with clip lock (the end with the two holes) [2a]
- 3) Introduction tube [3] with black end [3a]
- 4) Cleaning tube [4]
- 5) 0.1% Polyhexanide / 0.1% betaine cleaning spray (Prontosan[®] Wound Spray)

Fig. 1: Components of AlaxoLito Nasal Stent (left: AlaxoLito Nasal Stent, right: AlaxoLito Plus Nasal Stent)

3 Scope of Delivery

Contents of the AlaxoLito (Plus) Set:

- 1 AlaxoLito Box with:
 - 2 Nitinol braids
 - 1 Pusher
 - 1 Introduction tube
 - 1 Cleaning tube
 - 1 Bottle cleaning spray (Prontosan® Wound Spray 75 ml)
 - Instructions for use

Replacement Prontosan® Wound Spray bottles can be bought separately.

4 Use

The AlaxoLito Nasal Stent is a medical device and has been developed and tested for self-application by the patient. The nitinol braid consists of thin and sensitive metal wires so that the braid needs to be treated carefully and gently. The following instructions for handling and cleaning have to be obeyed in order to maintain the nitinol braids in good order over the allowed period of use of 18 months.

4.1 Preparation of the AlaxoLito stent for introduction

1. The nitinol braid is connected to the pusher [2] by fitting into each other the two parts [1a] and [2a] of the clip lock (see **Fig. 2**).
2. The other end of the pusher is inserted into the black end [3a] of the introduction tube [3] and the pusher is moved through the introduction tube. Then the nitinol braid is carefully and completely pulled into the introduction tube using the pusher whereby it self-compresses. The clip lock shall be located in the introduction tube. (see **Fig. 3**)

Fig. 2: Fitting into each other of the clip lock


Fig. 3: Insertion of the nitinol braid into the introduction tube

4.2 Introduction of the AlaxoLito stent into the nose

1. The introduction tube containing the nitinol braid is carefully introduced with the black end ahead into the left or right nostril as far as the clip lock still being outside of the nose (**Fig. 4**).
2. The pusher is held with the one hand and the introduction tube is slowly retracted over the pusher with the other hand. The nitinol braid self-expands and thereby attaches itself towards the walls of the nose. Thus the AlaxoLito stent is safely and stably positioned.
3. Once the introduction tube is fully removed the clip lock is pulled apart and thereby the pusher is separated from the nitinol braid.
4. Subsequently the other stent is introduced into the other nostril in the same way.
5. Cleaning of the introduction tube is described in *Section 4.4*.


Fig. 4: Scheme for positioning of the AlaxoLito Nasal Stent in the inferior nasal passage (left: AlaxoLito Nasal Stent, right: AlaxoLito Plus Nasal Stent)

4.3 Removal of the AlaxoLito stent

1. For removal the two nitinol braids are cautiously pulled out of the nose. Should the nasal mucosa have dried out it is recommended to first moisten the mucosa.
2. For that purpose a nasal spray containing sea salt, dexpanthenol and/or hyaluronic acid may be sprayed into the nose.

4.4 Cleaning of the AlaxoLito Nasal Stent

1. After retraction from the nose the **introduction tube** is rinsed in lukewarm (not hot !!!) flowing drinking water until all potentially adhering mucus is removed.

2. The outer surface of the introduction tube is sprayed with the cleaning spray (Prontosan® Wound Spray) which is also sprayed into the end of the introduction tube. The cleaning spray has to act for one minute and then the introduction tube is let dry.
3. The **nitinol braid** is rinsed in lukewarm (not hot !!!) flowing drinking water immediately after pulling out of the nose until all potentially adhering mucus is mostly removed.
4. The nitinol braid is pushed over the cleaning tube [4] and potentially adhering contaminations are carefully removed in lukewarm flowing drinking water with the fingers. For bringing the nitinol braid onto the cleaning tube its open end is cautiously pushed - obliquely fitted - over one end of the tube (do not pull!!!). Subsequently push the whole nitinol braid over the tube as shown in **Fig. 5**. Removal of contaminations from the nitinol braid with the fingers has to be done only in longitudinal direction, never in cross direction (i.e. not around the cleaning tube) as otherwise the sensitive nitinol braid may be damaged.


Fig. 5: Pushing of the nitinol braid over the cleaning tube

5. Once the nitinol braid is completely cleaned it is cautiously pulled off the cleaning tube.
6. Subsequently the nitinol braid is sprayed with the cleaning spray (Prontosan® Wound Spray) which has to act for one minute. Then the nitinol braid is rinsed once more in lukewarm flowing drinking water, again sprayed with the cleaning spray and let dry.
7. The **cleaning tube** is rinsed, too, in lukewarm flowing drinking water, sprayed with the cleaning spray which has to act for one minute and then let dry.
8. Should the **pusher** be contaminated it is also cleaned by rinsing in lukewarm flowing drinking water, spraying with the cleaning spray which has to act for one minute and let dry.

The three-step cleaning procedure for the AlaxoLito Nasal Stent serves for disinfection. By extensive rinsing the contaminations are efficiently removed. The 0.1% polyhexanide / 0.1% betaine cleaning spray (Prontosan® Wound Spray) additionally reaches all interspaces due to its modified surface tension. The ingredients possess excellent mucosa compatibility and a broad pathogen spectrum.

4.5 Storage of the AlaxoLito Nasal Stent

1. The dried nitinol braid shall not be stored in the introduction tube to prevent it from unnecessary pressure load.
2. After cleaning and drying all components of the AlaxoLito Nasal Stent are stored in the AlaxoLito Box until next use.
3. The AlaxoLito Set shall not be stored in direct sunlight and at temperatures below 4°C or above 50°C.

4.6 Period of use

1. The AlaxoLito Set after first use may be used for 18 months for nightly application.
2. After this period an increased risk of breakage of a nitinol wire due to material wear exists, especially for the round wires at the open end of the nitinol braid. Therefore, usage is no longer allowed.
3. The nitinol braid must be inspected for any damage prior to each use. If breakage of a wire is observed the nitinol braid must not be used anymore. A new AlaxoLito Set has to be used. Inspection is easiest performed on the cleaning tube during the cleaning procedure.
4. A cleaning spray bottle (Prontosan® Wound Spray) may be used for 12 months after first use. Expiry date see print on label.

5 Disposal

The AlaxoLito Set and the cleaning spray (Prontosan® Wound Spray) can be disposed of in domestic waste.

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