
# AlaxoStent

Gebrauchsanweisung  
Instructions for Use  
Manuel d'utilisation  
Istruzioni per l'uso  
Gebruiksaanwijzing

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## 1 Instructions

The AlaxoStent is a CE marked prescription medical device. It must only be used according to the briefing of the attending doctor.

### 1.1 CE Mark

The product bears the CE mark. The conformity evaluation assures that the product complies with the requirements of the EU directive for medical devices 93/42/EEG. According to the classification rules in Annex IX of the EU directive the AlaxoStent is a Class I product.

### 1.2 General Instructions

- The instructions for use is a component of the product. It has to be accessible to the user near to the product. Close adherence to these instructions is precondition for intended use and correct handling of the product.
- In case of damage resulting from use of third-party components or consumables no warranty claims apply.
- Alaxo is only liable for the product in respect of safety, reliability and functionality if it is used in compliance with the instructions for use.
- The instructions for use corresponds to the design of the product and the applicable safety norms as of going to print. All rights are reserved for the specified parts and components, processes, names and products.
- Full or partial reproduction subject to prior written approval by Alaxo.

### 1.3 General safety instructions

- The product has been developed and designed by Alaxo in such a manner that hazards during intended use are eliminated as far as possible. Nevertheless, we strongly recommend to strictly adhere to the following safety precautions in order to further reduce any residual risks to a minimum.
- It is not allowed to modify or reconstruct the product. In that case the CE mark lapses. Use of altered products is prohibited. In addition, Alaxo cannot take any warranty or liability for modified or reconstructed products. In order to assure safe application of the product users are liable for adherence to the instructions.
- The original packaging should be stored for any potential return shipment. Only the original packaging ensures optimal protection of the product during transport. Should return shipment become necessary during the warranty period Alaxo does not take any liability for damages resulting from insufficient packaging.
- The product must not be accessible to children.
- The user has to inspect the product for functional safety and good condition each time prior to application.
- The user has to be familiar with handling of the product.
- The cleaning instructions have to be observed and obeyed to. Other cleaning procedures must not be applied. The cleaning spray (Prontosan® Wound Spray) must only be used for disinfectant cleaning of the AlaxoStent and the accessories. Do not swallow!

### 1.4 Symbols and indications

- One Prontosan® Wound Spray bottle may be used for 12 months after first use.
- One nitinol braid may be used for 4 months after the first use.
- In the instructions for use, on packagings and on the product the following figures and symbols are used:

Consult Instructions for Use

Batch code

Month/Year

Use by: Month/Year

Order number

Manufacturer

CE Mark

Temperature limitation: store between 4° bis 50°C

Store inaccessible for children

## 1.5 Precautionary measures

The attending doctor has to explain handling of the AlaxoStent to the patient. Especially, the doctor has to explain correct application and positioning of the AlaxoStent in the throat and its removal.

The attending doctor has to advise on potential risks and complications connected to the AlaxoStent therapy and has to evaluate if the AlaxoStent therapy is applicable to the specific patient.

## 2 Product information

### 2.1 Intended use

- The AlaxoStent is a medical device and self-applied by the patient before going to sleep. The AlaxoStent is introduced via a nasal passage into the throat and positioned there in order to keep open the airway during sleep by mechanical splinting.
- The AlaxoStent is a prescription medical device. The attending doctor has to be sufficiently qualified in treatment of sleep apnea.
- Only one AlaxoStent is introduced through one nasal passage, never two or more at the same time.
- The AlaxoStent is clinically tested and CE marked for treatment of patients with obstructive sleep apnea or snoring.
- The attending doctor has to determine the disease state of the patient and if the AlaxoStent may be a useful therapy for this patient. The physician further has to examine if the nasal passages of the patient allow introduction of the AlaxoStent. The physician shall inform the patient which of the two nasal passages (left or right) is preferred for introduction of the AlaxoStent.
- In case of known lung diseases a pneumologic medical examination is required to evaluate if the AlaxoStent therapy can be applied to the specific patient.
- The AlaxoStent must be used only by one single patient.
- The AlaxoStent is suitable for repeated use up to 18 hours per application.
- The AlaxoStent must be used only by patients of age of 18 years and higher.
- The AlaxoStent is not suitable for life-supporting application in emergency care.

### 2.2 Contraindications

- In case of certain diseases the AlaxoStent should not be used or only with strong precautions. Any individual case requires thorough evaluation by the attending physician if the AlaxoStent therapy may be applicable.
- The AlaxoStent must not be used in case of inappropriate structure of the nasal passages, especially if leading to problems during insertion of the introduction tube into the nasal passage. Patients suffering from a hole in the nasal septum and similar anatomical disorders may not be accessible to treatment with the AlaxoStent.
- The AlaxoStent is not applicable for treatment of central apneas.
- The AlaxoStent shall not be used by persons with reduced mental ability, with reduced motor skills, with psychic restrictions, mental disorders or spastic syndromes.

### 2.3 Warnings

- Persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this stent.
- Persons with allergic reactions to polyhexanide or undecylenamidopropyl-betaine may suffer an allergic response to the cleaning spray.
- Do not use the AlaxoStent if this leads to pain! Do not terminate therapy without consulting your attending doctor.


### 2.4 Unintended use

- The only permitted use is splinting of the airway by introduction through a nasal passage into the throat. Other uses are prohibited.
- The AlaxoStent must not be introduced into orifices of the body other than the nose (e.g. ears or mouth).
- The AlaxoStent braid must not be pulled over any body parts or limbs.
- The product may only be used for the approved medical purposes. It must not be used for any other purposes, especially non-medical purposes.

## 2.5 Product description

The AlaxoStent Set comprises the following components (see **Fig. 1**):

- 1) Self-expanding nitinol braid [1] with clip lock [1a] (nitinol is a shape-memory metal, a nickel-titanium alloy widely used for medical devices)
- 2) Pusher [2] with clip lock (the end with the two holes) [2a]
- 3) Introduction tube [3]
- 4) Cleaning tube [4]
- 5) Fixation tape
- 6) 0.1% Polyhexanide / 0.1% betaine cleaning spray (Prontosan® Wound Spray)

**Fig. 1:** Components of AlaxoStent

## 3 Scope of Delivery

Contents of the AlaxoStent Start Set:

- 1 AlaxoStent Start-Box with:
  - 1 Nitinol braid (sealed in plastic tube)
  - 1 Pusher
  - 1 Introduction tube
  - 1 Cleaning tube
  - 1 Fixation tape
- 1 Bottle cleaning spray (Prontosan® Wound Spray 75 ml)
- Instructions for use

Contents of the AlaxoStent Annual Set:

- 1 AlaxoStent Annual Box mit:
  - 3 Nitinol braids (sealed in plastic tube)
  - 1 Pusher
  - 1 Introduction tube
  - 1 Cleaning tube
  - 1 Fixation tape
- 1 Bottle cleaning spray (Prontosan® Wound Spray 75 ml)
- Instructions for use


Replacement fixation tapes and Prontosan® Wound Spray bottles can be bought separately.

## 4 Use


The AlaxoStent is a medical device and has been developed and clinically tested for self-application by the patient. The nitinol braid consists of thin and sensitive metal wires so that the braid needs to be treated carefully and gently. The following instructions for handling and cleaning have to be obeyed to in order to maintain the nitinol braid in good order over the allowed period of use of 4 months.

#### 4.1 Preparation of the AlaxoStent for introduction

1. The plastic packaging containing the nitinol braid [1] is carefully opened with scissors at the end where the clip lock [1a] is positioned.
2. The nitinol braid is connected to the pusher [2] by fitting into each other the two parts [1a] and [2a] of the clip lock (see **Fig. 2**).
3. The other end of the pusher [2b] is inserted into black end [3a] of the introduction tube [3] and the pusher is moved through the introduction tube. Then the nitinol braid is carefully and completely pulled into the introduction tube using the pusher whereby it self-compresses. The clip lock shall be located in the introduction tube and the end of the stent shortly before the black end. (see **Fig. 3**)


**Fig. 2:** Fitting into each other of the clip lock

**Fig. 3:** Insertion of the nitinol braid into the introduction tube

#### 4.2 Introduction of the AlaxoStent into the throat

1. The introduction tube is carefully introduced with the black end ahead into the right or left nostril (tilt the head back!) and advanced horizontally through the nasal passage according to **Fig. 4** down to the throat. Exact positioning of the nitinol braid is explained to the patient by the prescribing doctor. In each case the AlaxoStent stays away from the epiglottis.
2. The smooth material and the bent form of the introduction tube serves for good steerability at the end of the nasal passage down to the throat and is optimally adapted to the human anatomy. Nevertheless, differences in simplicity of the introduction process may arise e.g. from anatomic variability from patient to patient.
3. In case of dry mucosal tissue in the nose pushing forward of the introduction tube may lead to increased friction. To wetten the nasal mucosa, after consultation of the attending physician, a nasal spray containing sea salt, dexpanthenol and/or hyaluronic acid may be applied horizontally into the nasal passage. The introduction tube smoothly moves forward along the wet mucosal tissue with low friction.

**Fig. 4:** Scheme for positioning of the AlaxoStent in the throat

4. After reaching the right position the pusher is held with the one hand and the introduction tube is slowly retracted over the pusher with the other hand. The nitinol braid self-expands first in the throat and then in the nasal passage and attaches itself towards the walls of the throat.
5. Once the introduction tube is fully removed the clip lock is pulled apart and thereby the pusher is separated from the nitinol braid. The AlaxoStent thus is in its position of use (see **Fig. 4**). Cleaning of the introduction tube is described in *Section 4.4*.
6. To secure the position of the AlaxoStent in the throat the part of the nitinol braid protruding from the nose is affixed with the fixation tape at the upper lip, the nose or the cheek. It has to be made sure that the fixation tape securely adheres to the nitinol braid.
7. Application of the AlaxoStent is completed and the patient can go to sleep.
8. Due to the opening force of the nitinol braid the throat is kept open and a collapse of the airway is prevented. Differences in the efficacy of splinting of the airway may be observed from patient to patient.

### 4.3 Removal of the AlaxoStent

After sleep the fixation tape is cautiously removed and the nitinol braid retracted from the nose and throat.

### 4.4 Cleaning of the AlaxoStent

1. After retraction from the nose the **introduction tube** is rinsed in lukewarm (not hot !!!) flowing drinking water until all potentially adhering mucus is removed.
2. The outer surface of the introduction tube is sprayed with the cleaning spray (Prontosan® Wound Spray) which is also sprayed into the end of the introduction tube. The cleaning spray has to act for one minute and then the introduction tube is let dry.
3. The **nitinol braid** is rinsed in lukewarm (not hot !!!) flowing drinking water immediately after retraction from the nose until all potentially adhering mucus is mostly removed.
4. The nitinol braid is pushed over the cleaning tube [4] and potentially adhering contaminations are carefully removed in lukewarm flowing drinking water with the fingers. For bringing the nitinol braid onto the cleaning tube its open end is cautiously pushed over one end of the tube (do not pull!!). Subsequently push the whole nitinol braid over the tube as shown in **Fig. 5**. Removal of contaminations from the nitinol braid with the fingers has to be done only in longitudinal direction, never in cross direction (i.e. not around the cleaning tube) as otherwise the sensitive nitinol braid may be damaged.

**Fig. 5:** Pushing of the nitinol braid over the cleaning tube

- Once the nitinol braid is completely cleaned it is cautiously pushed off the cleaning tube with the thumb and one finger without exerting any pressure from the lower end to the upper end (see **Fig. 6**). Please do not pull the closed end nor the clip lock of the nitinol braid as it will contract due to its design and cannot be pulled off the cleaning tube!

**Fig. 6:** Pushing off of the nitinol braid from the cleaning tube

- Subsequently the nitinol braid is sprayed with the cleaning spray (Prontosan<sup>®</sup> Wound Spray) which has to act for one minute. Then the nitinol braid is rinsed once more in lukewarm flowing drinking water, again sprayed with the cleaning spray and let dry.
  - The **cleaning tube** and the **pusher** are rinsed, too, in lukewarm flowing drinking water, sprayed with the cleaning spray (Prontosan<sup>®</sup> Wound Spray) which has to act for one minute and then let dry.
- The three-step cleaning procedure for the AlaxoStent serves for disinfection. By extensive rinsing the contaminations are efficiently removed. The 0.1% polyhexanide / 0.1% betaine cleaning spray (Prontosan<sup>®</sup> Wound Spray) additionally reaches all interspaces due to its modified surface tension. The ingredients possess excellent mucosa compatibility and a broad pathogen spectrum.

#### 4.5 Storage of the AlaxoStent

- The dried nitinol braid shall not be stored in the introduction tube to prevent it from unnecessary pressure load.
- After cleaning and drying all components of the AlaxoStent are stored in the AlaxoStent Box until next use.
- The AlaxoStent Set shall not be stored in direct sunlight and at temperatures below 4°C or above 50°C.

#### 4.6 Period of use

- The nitinol braid after first use may be used for four months for nightly application.
- After this period an increased risk of breakage of a nitinol wire due to material wear exists, especially for the round wires at the open end of the nitinol braid. Therefore, usage is no longer allowed.
- The nitinol braid must be inspected for any damage prior to each use. If breakage of a wire is observed the nitinol braid must not be used anymore. A new nitinol braid from the AlaxoStent Set has to be used. Inspection is easiest performed on the cleaning tube during the cleaning procedure.
- The pusher, the introduction tube, the cleaning tube and the cleaning spray (Prontosan<sup>®</sup> Wound Spray) may be used for 12 months after first use.
- After expiry of these allowed periods of use or if all nitinol braids contained in the AlaxoStent Set have been used up, a new AlaxoStent Set has to be employed.

## 5 Disposal

The AlaxoStent Set and the cleaning spray (Prontosan<sup>®</sup> Wound Spray) can be disposed of in domestic waste.

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