




COMPLETE PRODUCT CATALOGUE

Established in 1997. Proudly Canadian. 🇨🇦
Sold in over 80 countries worldwide.



Proudly Canadian and established in 1997, we specialize in the manufacturing of enzymatic solutions offered in different formats as well as the fabrication of silicone and plastic products. Our products are sold in over 80 countries worldwide and used primarily in surgery, endoscopy, and central sterilization areas.

Our in-house controlled environment production facility in Montreal, Canada adheres to strictly enforced quality and regulatory control systems. We are registered as a device manufacturing establishment with the United States FDA and we are in full compliance with CE and Health Canada Regulations. We continually invest in upgrading our plant, machinery, and processes to ensure high quality products at high volumes. Most products are available off-the-shelf.

At Batrik, we are committed to innovation and providing our customers with the finest quality products backed by world-class customer care.

SURGICAL PRODUCTS & GENERAL ACCESSORIES

04	SUTURE CLAMP TAGS
05	CLAMP JAW COVERS
06	VASCULAR SILICONE TIES
07	CAUTERY TIP CLEANER
08	SKINST [™] Surgical MULTI-TIP Markers
09	PIN GUARDS
10	SAFE-SOAK [™] SPAC & MATS
11	MAGNETIC MAT
12	TOURNY [™] + TOURNIQUETS
14	SHOOZCOVERS [™] SOFT LINE [™]
16	SHOOZCOVERS [™] TOUGH LINE [™]
18	SHOOZCOVERS [™] COATED LINE [™]
20	G.O.L.F.F. [™] ANTI-FOG SOLUTION
21	WTFOG [™] Reusable Anti-Fog Cloths
22	STACKABLE STEP STOOLS
23	BLUE CLOUD-GEL [™] ANTI-FATIGUE MATS
24	BLUE-CLOUD NIMBUS [™] ANTI-FATIGUE MATS

ENDOSCOPY PRODUCTS

25	G.O.L.F.F. [™] ANTI-FOG SOLUTION
26	BITE BLOCKS Endoscopy Bo2[™] / Bitestrong[™]
27	ENDOSCOPY BITESTRONG [™] BITE BLOCKS
28	COUVRETTE [™] VALVE COVERS
29	DAYBRUSH [™] CHANNEL BRUSHES
30	ARC [™] CHANNEL BRUSHES
33	VALVE / CONTROL BRUSHES
34	V-BRUSHES

INFECTION CONTROL & STERILE PROCESSING

35	TAG-IT [™]
39	SILICONE PIN MAT
40	CORNER POCKETS [™] & SOFT POCKETS [™]
41	LILA'S INSTRUMENT STAIN REMOVER
42	INSTRUMENT PROTECTORS
45	PrAmLi Power³ LUBE [™] INSTRUMENT LUBRICANT
46	S.T.I.C.S.

47	ID TAPE
48	TAPE-B-GONE [™] TAPE REMOVER
49	PADLOCKS SECURITY TAGS
53	SAVEFACE [™] FACE SHIELDS
54	SHIELD [™] FACE SHIELD
55	SCRUB PAD
56	Your First <input checked="" type="checkbox"/> [™] ENZYMATIC & NON-ENZYMATIC CLEANERS
57	Your First <input checked="" type="checkbox"/> [™] ENVIRO-SACHET [™]
58	Your First <input checked="" type="checkbox"/> [™] FOAMING SPRAY
59	Your First <input checked="" type="checkbox"/> [™] F1RST WIPES
60	Your First <input checked="" type="checkbox"/> [™] DRY SPONGES
61	Your First <input checked="" type="checkbox"/> [™] ENZYMATIC SPONGES & BEDSIDE SPONGE KITS
62	Your First <input checked="" type="checkbox"/> [™] NON-ENZYMATIC DETERGENT BEDSIDE SPONGE KITS
63	Your First <input checked="" type="checkbox"/> [™] TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

BRUSHES



64	ARC [™] WIRE STYLE BRUSHES
68	DOUBLE ARC [™] WIRE BRUSHES Asymmetrics / Symmetric
71	IN-BRUSH [™] ANTIMICROBIAL TWISTED WIRE BRUSHES
73	STAINLESS STEEL ARC [™] MICRO BRUSHES
74	INSTRUMENT CLEANING FLAT BRUSHES
75	ORTHOPEDIC RASPS CLEANING BRUSH
76	PRECISION BRUSH
77	INSTRUMENT CLEANING BRUSHES
78	ACETABULAR REAMER BRUSH
79	ALL PURPOSE REUSABLE BRUSHES
80	AUTOClave CLEANING BRUSH KITS

BATRIK SUTURE CLAMP TAGS

Non-toxic, radiopaque, and latex free, our Suture Clamp Tags are designed to protect and safely hold sutures without any damage.

Additionally, its foam base, affixed to a sterile area with the clamps, allows for easy insertion.

- Packed in a convenient foam base with easy peel-off adhesive backing, even while wearing gloves
- Available in sterile and non-sterile
- Available in 2, 3, or 5 pairs
- Pre-counted



OUR "EASY PEEL CORNER" IS A FEATURE UNIQUE TO BATRIK

PRODUCT SELECTION

ITEM # STERILE	ITEM# NON-STERILE	COLOR	DIMENSIONS	PACKAGING
11-1002	911-1002	Yellow	Standard ø (0.06" x 0.5" / 0.15 cm x 1.27 cm)	2 pairs / foam base
11-1003	911-1003	Yellow	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	3 pairs / foam base
11-1005	911-1005	Yellow	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1025	911-1025	Assorted	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1045	911-1045	White	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1065	911-1065	Blue	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1085	911-1085	Red	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-3002	911-3002	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	2 pairs / foam base
11-3003	911-3003	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	3 pairs / foam base
11-3005	911-3005	Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	5 pairs / foam base

STERILE PACKAGING: 1 Foam Base / Pouch – 5 Pouches / Box

NON-STERILE PACKAGING: Sold per foam base

BATRIK CLAMP JAW COVERS

Manufactured from medical grade silicone, our Clamp Jaw Covers reduce trauma to vessels requiring occlusion, particularly when a hemostatic clamp is utilized. The thread-like indents allow for gentle, constant and non-slip tension.

GOOD TO KNOW: OUR CLAMP JAW COVERS CONVERT ALMOST ANY INSTRUMENT TO ATRAUMATIC



- Gamma sterilized
- Available in 3 colors and 3 widths to fit most clamp sizes
- Radiopaque and latex free



PRODUCT SELECTION

ITEM # STERILE	ITEM # NON-STERILE	DIMENSIONS	COLOR
12-1000	912-1000	0.1" X 5" / 0.25 cm X 12.7 cm	Red
12-1001	912-1001	0.2" X 5" / 0.51 cm X 12.7 cm	White
12-1002	912-1002	0.3" X 5" / 0.76 cm X 12.7 cm	Blue

STERILE PACKAGING: 2 clamps per sterile pouch - 10 pouches per box

NON-STERILE PACKAGING: 20 clamps per bag (also available 100 units / bag)

BATRIK

VASCULAR SILICONE TIES

AVAILABLE IN 4 SIZES
AND 4 COLORS



Our Vascular Silicone Ties permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.


Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the BATRIK elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.

PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR
10-1000	0.075" X 16" (0.2 cm x 40 cm)	Mini-Red
10-1001	0.075" X 16" (0.2 cm x 40 cm)	Mini-Blue
10-1002	0.075" X 16" (0.2 cm x 40 cm)	Mini-Yellow
10-1003	0.075" X 16" (0.2 cm x 40 cm)	Mini-White
10-1010	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Red
10-1011	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Blue
10-1012	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Yellow
10-1013	0.120" x 16" (0.3 cm x 40 cm)	Maxi-White
10-1101	0.040" x 16" (0.1 cm x 40 cm)	Extra-Mini-Blue
10-1111	0.200" x 22" (0.5 cm x 56 cm)	Extra-Maxi-Blue
10-1113	0.200" x 22" (0.5 cm x 56 cm)	Extra-Maxi-White
10-3200	0.075" x 32" (0.2 cm x 81 cm)	Mini-Red
10-3201	0.075" x 32" (0.2 cm x 81 cm)	Mini-Blue
10-3202	0.075" x 32" (0.2 cm x 81 cm)	Mini-Yellow
10-3203	0.075" x 32" (0.2 cm x 81 cm)	Mini-White
10-3210	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Red
10-3211	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Blue
10-3212	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Yellow
10-3213	0.120" x 32" (0.3 cm x 81 cm)	Maxi-White



- Manufactured of medical grade silicone



- Do not twist on themselves when unpacked



- Easy to separate, eliminating the need the need for cardboard inserts



- Sterile, radiopaque, and latex free

PACKAGING: 2 ties per sterile pouch – 10 pouches per box

BATRIK CAUTERY TIP CLEANER

Our Cautery Tip Cleaner assists in removing material from cautery tips for both mono-polar and bi-polar electrodes.

- Foam pad with integrated abrasive material
- Adheres to disposable and non-disposable drapes
- Radiopaque
- Available exclusively in Non-Sterile
- Latex free



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
60-10000	Cautery Tip Cleaner Non-Sterile	1.8" x 1.8" / 4.5 cm x 4.5 cm	2,200 /Box

BATRIK SKINS™ SURGICAL MULTI-TIP SKIN MARKERS

Our SKINS™ Surgical Multi-Tip Skin Markers are designed with a beveled tip to enable the creation of fine or bold tracing lines by holding the pen at different angles.

- Available with or without 6 labels and a flexible 6" calibrated (in/mm) ruler
- Available in non-sterile for pre-operative use, or sterile for use in the operating room
- Surgical prep resistant
- Made of lasting gentian violet ink that is non-toxic and non-tattooing
- Latex free



PRODUCT SELECTION


ITEM #	TYPE	DESCRIPTION	PACKAGING
S-3000	Sterile	Multi-Tip fine marker with ruler and label set	50 / Box
S-3001	Sterile	Multi-Tip fine marker	50 / Box
S-30005	Sterile	Multi-Tip fine marker with ruler and label set (Volume Packaging)	500 / Box
S-30015	Sterile	Multi-Tip fine marker (Volume Packaging)	500 / Box
S-93000	Non-Sterile	Multi-Tip fine marker with ruler and label set	1000 / Box
S-93001	Non-Sterile	Multi-Tip fine marker	1000 / Box
S-94000	Non-Sterile	Ruler	100 / Bag
S-94001	Non-Sterile	Ruler (Bulk)	20,000 / Box
S-95000	Non-Sterile	Set of 3 labels	100 / Bag

SAVE by ordering our SKINS™ MARKERS in volume packaging. Cost effective and reduces waste.

BATRIK PIN GUARDS



- Safely cover the sharp ends of your K-Wire and Steinmann Pins
- Color coded by wire size which eliminates the risk of selecting of the wrong size
- A removable sleeve is provided to avoid injury and facilitate application (sterile only)
- Non-sterile model is available with or without plastic sleeves
- Packaging designed to reduce waste
- Latex free



PRODUCT SELECTION

ITEM # PIN GUARDS + SLEEVE	ITEM # PIN GUARDS + SLEEVE	ITEM # PIN GUARDS (NO SLEEVE)	COLOR	DIMENSIONS
14-1028	914-1028 - SLV	914-1028	Yellow	0.028" / 0.71 mm
14-1035	914-1035 - SLV	914-1035	Light Blue	0.035" / 0.89 mm
14-1045	914-1045 - SLV	914-1045	White	0.045" / 1.14 mm
14-1054	914-1054 - SLV	914-1054	Dark Red	0.054" / 1.40 mm
14-1062	914-1062 - SLV	914-1062	Green	0.062" / 1.57 mm
14-1078	914-1078 - SLV	914-1078	Pink	0.078" / 1.98 mm
14-1093	914-1093 - SLV	914-1093	Grey	0.093" / 2.36 mm
14-1109	914-1109 - SLV	914-1109	Beige	0.109" / 2.77 mm
14-1125	914-1125 - SLV	914-1125	Purple	0.125" / 3.17 mm
14-1141	914-1141 - SLV	914-1141	Brown	0.141" / 3.58 mm
14-1156	914-1156 - SLV	914-1156	Dark Blue	0.156" / 3.96 mm
14-1188	914-1188 - SLV	914-1188	Red	0.188" / 4.77 mm

STERILE PACKAGING: 80 units/box (2 units per sterile pouch - 40 pouches per box)

NON-STERILE PACKAGING: 100 units/box

BATRIK SAFE-SOAK™ SPAC

Our Safe-Soak™ SPAC floor suction disk quickly removes fluid waste from the operating room floor through the application of suction.



- Highly durable: it won't be damaged if stepped on
- Easily maneuverable by foot around operating room floor
- Latex free, disposable, and non-sterile
- The bottom opening permits fluid collection through a highly absorbent material that comes into contact with its surface
- The hose is securely affixed on its side to reduce risk of catching equipment
- Its highly absorbent base material allows for quick fluid wicking and removal when suction is applied

BATRIK SAFE-SOAK™ ABSORBENT FLOOR MATS

Our Safe-Soak™ Absorbent Floor Mats are designed to absorb a variety of fluids from the operating room floor quickly and efficiently, ensuring the floor beneath remains slip-free, safe, and dry during surgical procedures.

- Highly absorbent
- Facilitates quick and easy clean up
- Latex free, single-use, disposable, and non-sterile
- Conforms to uneven surfaces
- Skid-resistant backing (blue absorbent mat only)




PRODUCT SELECTION

ITEM #	DESCRIPTION	ABSORBENCY	SIZE	PACKAGING
SS-20-1000	Safe-Soak™ / White extra absorbent mat	4.5 L / mat	28" x 48" (71 cm x 122 cm)	20 / Box
SS-20-2000	Safe-Soak™ / Blue absorbent mat with skid-resistant backer	2.5 L / mat	28" x 40" (71 cm x 102 cm)	20 / Box
SS-20-3000	Safe-Soak™ / Yellow absorbent mat	3 L / mat	30" x 52" (76 cm x 132 cm)	15 / Box
SS-21-1000	Safe-Soak™ SPAC	N/A	8.625" (22 cm diameter)	12 / Box

BATRIK MAGNETIC MAT

Our Magnetic Mat provides a sterile, neutral field for the safe placement and hands-free transfer of instruments. It is latex free, reusable, and autoclavable (up to 134° C).



PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
MAG-02-1	Magnetic Mat	12" x 15" / 30 cm x 38 cm	1 Unit
MAG-02-5	Magnetic Mat (Volume Packaging)	12" x 15" / 30 cm x 38 cm	5 Units

*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy for reprocessing of medical devices in a healthcare setting

BATRIK

TOURNY™ + TOURNIQUETS

Our Tourniquets are comfortable and easy to use on all patients, including the elderly, undernourished, children, and newborns.



AVAILABLE COLORS



- Available in different sizes including compact sizes, which are convenient for phlebotomy trays
- Innovative latex free material has a more comfortable texture for patients, puts less stress on skin, and safeguards hospital staff and patients from latex allergens

SINGLE-USE BENEFITS

- Reduced risk of cross-contamination (patient to patient or patient to healthcare provider)
- Reduced risk of spread of nosocomial infections
- May improve your facility's compliance with OSHA's 1910-1030 (d)(1) "Blood-borne Pathogens Standard"



FOR YOUR CONVENIENCE:
 25 UNITS PER ROLL
 COME IN A DISPENSER BOX
 FOR EASY ACCESS



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	DIMENSIONS	PACKAGING
K-002	TOURNY™ + Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 / Bag
K-004	TOURNY™ + Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	100 / Box
K-009	TOURNY™ + Non-Sterile	Blue	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
K-010	TOURNY™ + Non-Sterile	White	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
K-400	TOURNY™ + Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
K-400PB	TOURNY™ + Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
K-100	TOURNY™ + Non-Sterile	Assorted (3 rolls of each) Blue, Orange, Green, and Pink	1" x 18" / 2.5 cm x 45 cm	12 rolls 300 / Box
K-254	TOURNY™ + Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	250 / Box
K-2-18RL	TOURNY™ + Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-2-18PBRL	TOURNY™ + Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-3-18RL	TOURNY™ + Non-Sterile	Green	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-6-18RL	TOURNY™ + Non-Sterile	Pink	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-7-18RL	TOURNY™ + Non-Sterile	Orange	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-2-29RL	Long TOURNY™ + Non-Sterile	Blue	1" x 29.5" / 2.5 cm x 75 cm	1 Roll

SHOOZCOVERS™ | SOFT LINE™


Batrik ShoozCovers™ Soft Line™ is manufactured in Canada from high quality 35 g/m² polypropylene fabric and 100% latex free elastics for performance and durability you can rely on. Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.

SHOOZ-REPELLENT™ TECHNOLOGY:

The **Soft Line™** integrates Shooz-Repellent™ technology which creates a protective barrier that repels fluids, delivering basic incidental fluid contact repellent. Shoe covers are not waterproof or water-resistant.

ANTI-SKID TREAD:

Our anti-skid tread reduces the risk of slipping with **triple ultrasonic welds** that increases tear strength and delivers easy on-and-off convenience.



MADE IN CANADA 



OUR SOFT LINE™ IS PERFECT FOR:



MEDICAL
ESTABLISHMENTS



PHARMACEUTICAL
LABS



GMP
SITES



ELECTRONIC
MANUFACTURING

LOLLY BLUE | **SOFT LINE™**



WHY YOU'LL LOVE THE SOFT LINE™:

- Both shoe covers are manufactured in Canada with top-quality 35 g/m² polypropylene fabric and 100% latex-free USA-made elastics¹.
- This line is particle count tested and complies with IEST-RP-CC003.4.
- Other colors available upon request CC003.4.



AVAILABLE IN TWO VERSATILE SIZES

BUBBLE GUM PINK | **SOFT LINE™**



SUPPORT A CAUSE, WEAR PINK:

- A portion of our Bubble Gum Pink shoe cover sales will be donated to charity to help further breast cancer research. Join us in making a meaningful contribution by choosing Pink ShoozCovers™!



AVAILABLE IN UNIVERSAL

PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING
SC-1210-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 1,000 (10 bags of 100 units)
SC-1216-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)
SC-1226-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 600 (6 bags of 100 units)
SC-1220-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 1,000 (10 bags of 100 units)
SC-1616-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Bubble Gum Pink	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags

¹Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.

SHOOZCOVERS™ | TOUGH LINE™

EXTRA-LARGE FOR EXTRA-TOUGH WORK:


Our Batrik ShoozCovers™ Tough Line™ comes in an **extra-large** size that fits over a **men's size 15 steel-toed boot**, delivering optimal full coverage and protection while you work.¹

ENHANCED DURABILITY:

These shoe covers are formulated with our **thickest fabric yet!** Made of **~55 g/m² polypropylene fabric**, these highly durable covers provide added protection, increase tear strength and reduce linting.



SHOE SIZES 12-15 (45-49.5 EU)
~ 47 CM X 19 CM / 18.5" X 7.5"



MADE IN CANADA 



OUR TOUGH LINE™ IS PERFECT FOR:



HOMES OR
COMMERCIAL
BUILDINGS



SHIELDING
AGAINST MUD
& GREASE



MEDICAL
ESTABLISHMENTS




PHARMACEUTICAL
LABS

¹ SHOOZCOVERS™ ARE MADE LONGER (42 & 47CM) FOR ACCURATE SHOE SIZE SELECTION.

BLACK LIQUORICE | TOUGH LINE™

WHY YOU'LL LOVE THE TOUGH LINE™:

- Triple ultrasonic welds that are long lasting to reduce waste
- High quality made in the USA latex-free elastics are triple reinforced to contract the shoe cover around your foot and deliver a strong and secure fit²
- Shooz-Repellent™ technology delivers basic incidental fluid contact repellency³
- Other colors available upon request


XL

AVAILABLE IN EXTRA LARGE

ANY TIME, ANYWHERE:



RESIDENTIAL



COMMERCIAL



INDUSTRIAL

PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING
SC-1922-PB	ShoozCovers™ - Tough Line™ ~55 g/m ² polypropylene non-woven fabric	Black Liquorice	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 200 (2 bags of 100 units)

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags

² Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product. ³ Shooz-Repellent™ technology provides basic incidental fluid contact repellent. It is not waterproof nor water resistant.


SHOOZCOVERS™ | COATED LINE™

STAY DRY:

Our Batrik ShoozCovers™ Coated Line™ is made from a ~55 g/m² blend of **polypropylene** and **polyethylene**, providing **high-quality water resistance**.¹ Ideal for environments where **maximum protection against fluids** is required.

TRIPLE WELDED TECHNOLOGY:

These shoe covers are reinforced with **triple ultrasonic welds**, a technique that joins materials together without the need for adhesives or stitching. This technique **increase tear strength and reduce linting**, for long lasting, resilient covers.



MADE IN CANADA 



OUR COATED LINE™ IS PERFECT FOR:



MEDICAL
ESTABLISHMENTS



PHARMACEUTICAL
LABS



GMP SITES




ELECTRONIC
MANUFACTURING

¹OUR SHOOZCOVERS™ COATED LINE™ IS RESISTANT AGAINST ALL NON PETROLEUM BASED LIQUIDS.

MARSHMALLOW WHITE | COATED LINE™

WHY YOU'LL LOVE THE COATED LINE™:

- **Superior protection using our unique blend of fabrics** which protects against not only water but all non-petroleum-based liquids, making them the ideal choice for environments where spills and leaks are a concern
- **Elastics around the ankle create a strong and secure fit** which is especially important when working in hazardous environments where a misplaced shoe cover could lead to contamination or injury
- **Smart design** with a large opening, providing a better fit and easy on-and-off convenience
- **Other colors available upon request**



UN XL

AVAILABLE IN TWO VERSATILE SIZES

PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING
SC-2212-NPB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)
SC-2212-NPDB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 boxes of 150 units)
SC-2122-NPB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)
SC-2122-NPDB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 boxes of 120 units)

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags

BATRIK

WTFOG™ by G.O.L.F.F.™

REUSABLE ANTI-FOG CLOTHS

"WHAT THE FOG?!"

We've all experienced it.

Eyewear fogging up while wearing a mask is not only unpleasant and highly frustrating, but it can be dangerous, too.

INTRODUCING THE FOG FIX:

Our WTFog™ reusable anti-fog cloths prevent fogging on all types of eyewear. Ideal for use when wearing a mask or face shield, WTFog™ takes you from "What The Fog?!" to confidently ready for your next patient in a snap.

WHY YOU'LL LOVE IT:

- Eliminates fogging on all glass, polycarbonate, and plastic eyewear lenses
- Each reusable nano-infused cloth lasts up to 600 wipes or 6 months
- Microfiber cloths are infused with our proprietary heavy duty anti-fog treatment
- Quick and easy application - just wipe and wear
- Lint-free, scratch-free, no greasy smearing
- Long lasting protective coating provides effective protection against fog, humidity, sweat, and fast temperature changes.
- 100% Natural Titanium Dioxide (TiO2) Nano Formula developed in Germany meets strict REACH standards
- Non-Sterile



"To be able to complete a whole appointment without adjusting my eye protection makes maintaining COVID PPE protocols that much easier."

DR. MOLLY FRIEDMAN

HOLLAND BLOORVIEW KIDS REHABILITATION HOSPITAL
TORONTO, ON

PERFECT WHEN:

- You wear glasses/goggles and a mask/face shield
- Your profession regularly exposes you to fogging and you need crystal clear vision without the risk or hassle of fogging up



DIRECTIONS FOR USE:

Open resealable pouch and remove the WTFog™ reusable anti-fog cloth. Wipe cloth on clean surface (glasses/goggle lens) being careful not to remove the protective coating delivered by the cloth. There is no wait time required before wearing your protected item. Enjoy your fog-free experience!

CLOTH STORAGE & CARE:

After use, keep cloth in the resealable pouch until next use, ensuring the seal is fully closed. If cloth becomes too dry, add up to 1 ml of water to pouch to reactivate cloth. For best results, store at ambient temperature/humidity. Dispose cloth after 600 uses or 6 months.

External surface use only.
Keep out of reach of children.
Do not wash

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
WTFOG-30	30 WTFog™ in a dispenser box measuring: 3 7/16" x 5 2/16" x 6 2/16" (8.7 cm x 13 cm x 15.4 cm)	30 WTFog™ / dispenser box
WTFOG-100	100 WTFog™ - bulk packed	100 WTFog™ / box
WTFOG-1200	1200 WTFog™ - bulk packed	1200 WTFog™ / case (12 boxes of 100 bulk packed)

1 WTFog™ Cloth per Resealable Pouch

WTFog™ Cloth: 5.5" x 5.5" (13.9 cm x 13.9 cm)

WTFog™ Resealable Pouch: 3.5" x 5" (8.9 cm x 12.7 cm)

IMPERIAL SURGICAL LTD.

STACKABLE STEP STOOLS

Made of lightweight cast aluminium and protected with an antimicrobial powder coating, Imperial Surgical Ltd.'s patented Stackable Step Stools have been designed with infection prevention as a prime focus.



PATENTED

- Weight of stool: 8.3 lbs (3.8 kg)
- Dimensions: 17.4" L x 13" W x 4.7" H (44.2 cm L x 33 cm W x 11.9 cm H)
- Antimicrobial coating prevents the spread of dangerous micro-organisms and suppresses the growth of mold, bacteria, yeast, and mildew
- Cart washable. Model OR-36364 (without mat) can be sterilized
- Anti-slip surface for enhanced stability
- Rubber-free corners avoids fluid collection
- No welds or large hollow openings to trap debris or fluids
- Transport cart available
- 300 lb (136 kg) weight capacity
- 7 year warranty covers all parts (step stool only)

PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR
OR-36363	Stackable Step Stool with mat	Black
OR-36363-02	Stackable Step Stool with mat	Dark Blue
OR-36363-03	Stackable Step Stool with mat	Green
OR-36363-07	Stackable Step Stool with mat	Orange
OR-36364	Stackable Step Stool without mat	N/A

PACKAGING: Sold by Unit

OPTIONAL ACCESSORIES

ITEM #	DESCRIPTION	PACKAGING
OR-3693	Transport cart	1 Unit
809-1705	Reusable Anti-Fatigue Mat with adhesive strips to secure it in place	5 / Box
899-1710	Single-Use Anti-Fatigue Mat with adhesive strips to secure it in place	10 / Box

REPLACEMENT PARTS

ITEM #	DESCRIPTION	PACKAGING
OR-36363-1A	Replacement mat for OR-36363	1 Unit
OR-36363-1B	Replacement Foot Pegs	1 Unit

TIP: COLOR-CODE YOUR STEP STOOLS BY DEPARTMENT



BATRIK BLUE CLOUD-GEL™ ANTI-FATIGUE MATS


Our Blue Cloud-Gel™ Anti-Fatigue Mats have been scientifically engineered for use in the O.R. using patented gel technology and premium foam to provide the ultimate comfort and support for surgeons, nurses, and technicians particularly when movement is restricted for prolonged durations.

- Proprietary 3/4" thickness Blue Cloud-Gel™ foam technology will not "bottom-out" like conventional foam mats
- Conforms to and supports feet promoting efficient circulation, reducing physical discomfort and fatigue
- Promotes proper posture while diminishing back, knee, and foot discomfort
- Latex free, water resistant, and nonabsorbent
- Reusable, easy to clean, and Cart Washer compatible (contact us for more details)
- Excellent friction properties when damp or dry surgical booties are worn
- Completely sealed with an RF weld for excellent protection against moisture absorption
- US Standards: CALL 117 and NFPA 260 Compliant
- 1 Year Warranty

DID YOU KNOW? STUDIES HAVE DEMONSTRATED THAT ANTI-FATIGUE MATS IMPROVE PERFORMANCE, ENERGY, AND POST-OPERATIVE DISCOMFORT.

PRODUCT SELECTION


ITEM #	DIMENSIONS	PACKAGING
800-2401	18" X 24" / 46 cm x 61 cm	1 / Box
800-2405	18" X 24" / 46 cm x 61 cm	5 / Box
800-3201	20" X 32" / 51 cm x 81 cm	1 / Box
800-3203	20" X 32" / 51 cm x 81 cm	3 / Box
800-4801	20" X 48" / 51 cm x 122 cm	1 / Box
800-4802	20" X 48" / 51 cm x 122 cm	2 / Box
800-7201	30" X 72" / 76 cm x 183 cm (Special Order)	1 / Box



18" X 24"



20" X 32"




20" X 48"

BATRIK BLUE-CLOUD NIMBUS™ ANTI-FATIGUE MATS

Our Blue-Cloud Nimbus™ Anti-Fatigue Mats are designed for general hospital use. Made with proprietary high-density “energy-return” polyurethane materials, they provide unmatched comfort, relief, and ergonomic support.

- Our Energy-Return properties will help reduce leg, foot, lower back discomfort, and work-related injuries
- Longer lasting support as 3/4” thickness proprietary foam blend will not bottom out like most other foam based anti-fatigue mats
- Antimicrobial properties blended into the mat provides added protection
- Highly water resistant surface will not absorb water or dirt
- Lightweight and easy to clean
- Textured design on surface and high-traction base ensures greater safety. (US Standard NFSI 101-C)
- Beveled edging reduces tripping hazards and anti-slip backing enhances safety
- Durable design is ideal for commercial applications, yet still ergonomically and aesthetically pleasing for home or office environments
- Flammability/Burn resistant per CFR 16 Parts 1630/1631
- Made with eco-friendly materials, is non-toxic, and latex free
- 5 Year Warranty (non-industrial usage), 4 Year Warranty (industrial usage)

GOOD TO KNOW: IF WITHIN
30 DAYS YOU ARE NOT 100%
SATISFIED, CONTACT YOUR
DEALER FOR A REFUND
(EXCLUDES TRANSPORT,
HANDLING, AND TAXES)



Available in
Taupe or Black

PRODUCT SELECTION

ITEM #	DIMENSIONS INCHES / CM	PACKAGING
BLACK COLOR		
809-3201	20" X 32" / 51 cm x 81 cm	1 / Box
809-3204	20" X 32" / 51 cm x 81 cm	4 / Box
809-4801	20" X 48" / 51 cm x 122 cm	1 / Box
809-4804	20" X 48" / 51 cm x 122 cm	4 / Box
809-3601	24" X 36" / 61 cm x 91 cm	1 / Box
809-7201	20" X 72" / 51 cm x 183 cm	1 / Box
809-6001	36" X 60" / 91 cm x 152 cm	1 / Box
TAUPE COLOR		
810-3201	20" X 32" / 51 cm x 81 cm	1 / Box
810-3204	20" X 32" / 51 cm x 81 cm	4 / Box
810-4801	20" X 48" / 51 cm x 122 cm	1 / Box
810-4804	20" X 48" / 51 cm x 122 cm	4 / Box
810-3601	24" X 36" / 61 cm x 91 cm	1 / Box
810-7201	20" X 72" / 51 cm x 183 cm	1 / Box
810-6001	36" X 60" / 91 cm x 152 cm	1 / Box

BATRIK
G.O.L.F.F.™ ANTI-FOG SOLUTION

Our G.O.L.F.F.™ Anti-Fog Solution solves the fogging problem often encountered during Laparoscopic surgical procedures.

- Prevents the occurrence of fog on lenses
- Ensures clearer vision and avoids unnecessary interruptions during surgery
- Set includes an easy-squeeze 6 mL bottle and an adhesive radiopaque sponge for secure placement
- OEM (available)
- Sterile
- Latex free

**PRODUCT SELECTION**

ITEM #	DESCRIPTION	PACKAGING
G.O.L.F.F.	Anti-Fog Solution & Adhesive Sponge set	20 / Box

BATRIK

ENDOSCOPY BO2™ BITE BLOCKS


WITH NASAL OXYGEN DELIVERY

Our BO2™ Bite Blocks significantly enhance patient safety by providing supplemental oxygenation during gastrointestinal endoscopy, preventing hypoxia in sedated patients, and ensuring superior patient comfort.

- Saves time and money by eliminating the need for separate nasal prongs and oxygen masks
- Designed to be free of sharp edges to ensure patient comfort and safety
- Available exclusively in Single-Use
- Optional retention strap is available to ensure a secure fit
- The Bite Blocks and optional retention straps are all latex free, lead free and DEHP free

BO2
(FRONT SIDE)

BO2
(BACK SIDE)



PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
BO2	Single-Use	18mm / 54FR, includes strap	100 / Box
BO2-N	Single-Use	18mm / 54FR, without strap	100 / Box

BATRIK ENDOSCOPY BITESTRONG™ BITE BLOCKS

ADULT & PEDIATRIC BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.

Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

SINGLE USE

- Available without dental retention

REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap

PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
B1	Single-Use Pediatric	16 mm / 48FR, includes strap	100 / Box
B1-25	Single-Use Pediatric	16 mm / 48FR, includes strap	25 / Box
B1-N	Single-Use Pediatric	16 mm / 48FR, without strap	100 / Box
B3	Single-Use	20 mm / 60FR, includes strap	100 / Box
B3-N	Single-Use	20 mm / 60FR, without strap	100 / Box
B3-SA	Single-Use	20 mm / 60FR, with strap <u>attached</u> to one side	100 / Box
B3-NDR	Single-Use	20 mm / 60FR, no dental retention rim	100 / Box
B3-NDR-N	Single-Use	Same as B3-N, no dental retention rim, without strap.	100 / Box
B4	Single-Use	19 mm / 60FR	100 / Box
B4T	Single-Use	19 mm / 60FR, includes strap	100 / Box
B4T-N	Single-Use	19 mm / 60FR, without strap	100 / Box
B6	Reusable	19 mm / 60FR	100 / Box
B6T	Reusable	19 mm / 60FR, includes strap	100 / Box
B6T-N	Reusable	19 mm / 60FR, without strap	100 / Box
B6-25	Reusable	19 mm / 60FR	25 / Box
B6T-25	Reusable	19 mm / 60FR, includes strap	25 / Box
B6T-N-25	Reusable	19 mm / 60FR, without strap	25 / Box
BS	Single-Use	Bite-Block straps 21.7" x 0.79" / 55 cm x 2 cm	100 / Box



BS



B1



B3



B4



B6

BATRIK

COUVRETTE™ VALVE COVERS

Our Couvrette™ Valve Covers eliminate the need for manual cleaning and reprocessing of reusable valves. Our single use Valve Covers effectively reduce cross contamination between patients and also reduce fluid leakage from the parts during procedures.

- Easily and securely fit to all valves
- Reduce the potential for errors
- Single-Use
- Non-Sterile
- Latex free



16-0002

16-0004

TIP: REUSABLE VALVE COVERS ARE NOTORIOUSLY DIFFICULT TO CLEAN.

SAVE TIME AND REDUCE RISK WITH COUVRETTE™.

COMING SOON!



16-0022

16-0044

PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	PACKAGING
16-0002	Biopsy Valve Cover - Olympus®	Blue	100 / Box
16-0004	Biopsy Valve Cover - Pentax®	Red	100 / Box
16-0022	Biopsy Valve Cover - Olympus® & Fujinon®	Light Blue	100 / Box COMING SOON!
16-0044	Biopsy Valve Cover - Pentax®	Pink	100 / Box COMING SOON!


NOTE: Contact us to verify compatibility with your scope

BATRIK DAYBRUSH™ CHANNEL BRUSHES

REUSABLE

Our Reusable DayBrush™ Channel Brushes can be used for cleaning endoscopes for an entire day. In fact, it's been proven to effectively clean up to 10 scopes before the need to discard.*

- Soft tip ensures easy manipulation within the channel with no risk of damage
- The back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse.
- Can be steam sterilized.
- Bristle volume is 20mm in length – nearly double most brushes on the market. A large brush provides more surface contact during the cleaning process.
- Significant per use savings when compared to other reusable or disposable channel cleaning brushes.



() Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications, the brush may have the appearance that it can continue to clean effectively, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.*

PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM

ITEM #	COLOUR	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
30-18050	White	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
30-18050-T	Yellow	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
30-18050-2T	Yellow	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double
30-23050	White	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
30-23050-T	Yellow	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
30-23050-2T	Yellow	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double

PACKAGING : 50 / BOX

OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.


"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

BATRIK ARC™ CHANNEL BRUSHES SINGLE-USE

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, materials, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, simply call us. We have over 50 different brushes available and can custom fabricate to your precise specifications.



OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

SERIES 311 & 333



PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
311-10050	39" / 100 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
311-18050	70" / 180 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
311-23050	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
311-23050-2	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Double

PACKAGING: 50 / BOX

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
333-23056	90" / 230 cm	0.94" / 24 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Single
333-23056-2	90" / 230 cm	0.94" / 24 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Double

PACKAGING: 50 / Box

SERIES 300 & 322



PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
300-18050	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
300-23050	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
300-23050-2	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double


PACKAGING: 50 / Box

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
322-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.07" / 1.75 mm	Single
322-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.07" / 1.75 mm	Double

PACKAGING: 50 / Box

SERIES 344




PRODUCT SELECTION

3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM - METALLIC PURPLE DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
344-10030	39" / 100 cm	0.79" / 20 mm	0.12" / 3 mm	0.06" / 1.55 mm	Single
344-18030	70" / 180 cm	0.79" / 20 mm	0.12" / 3 mm	0.06" / 1.55 mm	Single
344-23030	90" / 230 cm	0.79" / 20 mm	0.12" / 3 mm	0.06" / 1.55 mm	Single
PACKAGING: 50 / Box					

SERIES 35



PRODUCT SELECTION

2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
35-12025	47.2" / 120 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single
35-23025	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single
35-23025-2	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Double
PACKAGING: 50 / Box					

2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
35-12020	47.2" / 120 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
35-18020	70.8" / 180 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
35-23020	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
35-23020-2	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Double
PACKAGING: 50 / Box					

BATRIK

VALVE / CONTROL BRUSHES

27-05197



5.00" / 13 cm

27-57787



6.00" / 15 cm

40-18511



8.07" / 20.5 cm

36-10010



4.3" / 11 cm

27-15311




4" / 10.3 cm

PRODUCT SELECTION

ITEM #	DESCRIPTION	LENGTH	BRUSH 1 LENGTH	BRUSH 2 LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)
27-05197	Double Trumpet 5/7 White Color	5.00" / 13 cm	0.20" / 5 mm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-05197-50						
27-57787	Double. Universal 5/11 Yellow Color	6.00" / 15 cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-57787-50						
40-18511	Double Universal 5/11 Extra Long Blue Color	8.07" / 20.5 cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (reusable and can be steam sterilized). Volume Packaging: Bag of 50 use ITEM # 40-18511-50						
36-10010	Single end large Red Color	4.3" / 11 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 36-10010-50						
27-15311	Single end large Blue Color	4" / 10.3 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-15311-50						

BATRIK V-BRUSHES EXTERNAL CLEANING BRUSHES

Our External Cleaning Brushes for endoscopes / optical instruments provide excellent cleaning capabilities. With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments. Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.



PRODUCT SELECTION

ITEM #	DESCRIPTION	BRUSH LENGTH / (Ø)	HANDLE LENGTH / (Ø)	TOTAL LENGTH	PACKAGING
36-20005	Single-Use External Cleaning Brush	70 mm / 25 mm	80 mm / 3 mm	150 mm	50 / Box
36-20010	Reusable External Cleaning Brush	90 mm / 25 mm	120 mm / 12 mm	210 mm	3 / Bag

BATRIK

TAG-IT™

Batrik Tag-It™ is our new line of information tags developed to help medical professionals identify and document a wide range of instructions. Whether you need to identify instruments, trays or carts, both before or after sterilization, Tag-It™ has got you covered. These durable and tear-resistant tags are perfect for indicating special handling instructions.

WHY YOU'LL LOVE THEM:

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C) and are compatible with Gas plasma sterilization (STERRAD®100 NX Sterilization System).
- The synthetic paper is durable and will not wrinkle or absorb water. Made from recyclable plastic promoting environmental sustainability.
- Robust and easy to use tie wraps designed to securely attach to a wide range of trays, instruments, and other medical devices.
- Available in over 8 languages.

KEY PRODUCT FEATURES:

- **High Quality Permanent Marker**
- **Steam or Plasma indicator Dots** provide confirmation of processing
- **Eco-friendly, Waterproof & Tearproof Tags**
- **Large Quantities:** Each bundle of Tag-It™ comes with 100 versatile tags, to quickly identify instruments that require attention.
- **Available in multiple linguistic choices**

**TI-3503**

All tags available with indicator dots
Provides visual confirmation that the tag was processed

**PRODUCT SELECTION:**

ITEM #	COLOR	DESCRIPTION	SIZE	PACKAGING
TI-3502	Blue	VENDOR	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3503	Green	LOANER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3504	Red	PRIORITY	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3505	Yellow	RE-RUN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3506	Pink	MISSING	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3507	Orange	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3509	Black	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B1	White	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B4	Red	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35C22	Light Blue	VENDOR CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35L3	Green	LOANER CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35T2	Blue	QUICK TURNOVER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-PM09	Black	FINE PERMANENT MARKERS	1.5" Ø x 5.5" / 3.81cm Ø x 12.7 cm	36 / Bag

All information tags are single-use. **PACKAGING:** 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for steam dots.

PROCESSING INDICATOR DOTS: Add suffix "PL" to end of code for plasma dots.

PRODUCT SELECTION:

VENDOR / BLUE

ITEM #	SIZE
TI-3502	3" x 5" (8 cm x 13cm)



VENDOR

Company Name _____

Tray Name _____

Date of Surgery _____

WWW.BATRIK.COM TAG-IT™ TI-3502

BLANK / WHITE


ITEM #	SIZE
TI-35B1	3" x 5" (8 cm x 13cm)



WWW.BATRIK.COM TAG-IT™ TI-35B1

LOANER / GREEN

ITEM #	SIZE
TI-3503	3" x 5" (8 cm x 13cm)



LOANER

Loaned From _____

Loaned To _____

By _____ Date _____

WWW.BATRIK.COM TAG-IT™ TI-3503

BROKEN / RED

ITEM #	SIZE
TI-35B4	3" x 5" (8 cm x 13cm)



BROKEN DO NOT USE

DATE: _____ PHONE: _____

DEPARTMENT/CONTACT: _____

ITEM: _____

CRITICAL: (必) Yes No

INSTRUMENT SET: _____

ASSEMBLED BY: _____

WWW.BATRIK.COM TAG-IT™ TI-35B4

PRIORITY / RED

ITEM #	SIZE
TI-3504	3" x 5" (8 cm x 13cm)



PRIORITY

Room # _____

Time Needed _____

Initials _____

WWW.BATRIK.COM TAG-IT™ TI-3504

VENDOR IN + OUT / LIGHT BLUE

ITEM #	SIZE
TI-35C22	5" x 3" (13 cm x 8 cm)



VENDOR

Check In

Tray Received From: _____

Tray Name: _____

Received By: _____

Date: _____ Surgery Date: _____

☐ CPD Notified Prior to Receipt of Trays

☐ Received >24 Hours Before Case

☐ Inventory List Available

☐ Inventory Check Complete

☐ All Instruments in Good Condition

☐ Containers in Good Condition

☐ All Trays Labeled (Patient & Surgeon)

☐ All Trays <25 pounds

☐ Written Instructions for Use Available

☐ Disassembled/Decontaminated

WWW.BATRIK.COM TAG-IT™ TI-35C22

VENDOR

Check Out

☐ Disassembled/Decontaminated

☐ Inventory Check Complete

☐ All Items Complete and Accounted For

Returned By: _____

Received By: _____

Date: _____

Notes: _____

WWW.BATRIK.COM TAG-IT™ TI-35C22

RE-RUN / YELLOW

ITEM #	SIZE
TI-3505	3" x 5" (8 cm x 13cm)



RE-RUN

Room # _____

Time Needed _____

Initials _____

WWW.BATRIK.COM TAG-IT™ TI-3505

MISSING / PINK

ITEM #	SIZE
TI-3506	3" x 5" (8 cm x 13cm)



MISSING

Room # _____

Time Needed _____

Initials _____

WWW.BATRIK.COM TAG-IT™ TI-3506

LOANER IN + OUT / GREEN

ITEM #	SIZE
TI-35L3	5" x 3" (13 cm x 8 cm)



LOANER

Check In

Tray Received From: _____

Tray Name: _____

Received By: _____

Date: _____ Surgery Date: _____

☐ CPD Notified Prior to Receipt of Trays

☐ Received >24 Hours Before Case

☐ Inventory List Available

☐ Inventory Check Complete

☐ All Instruments in Good Condition

☐ Containers in Good Condition

☐ All Trays Labeled (Patient & Surgeon)

☐ All Trays <25 pounds

☐ Written Instructions for Use Available

☐ Disassembled/Decontaminated

WWW.BATRIK.COM TAG-IT™ TI-35L3

LOANER

Check Out

☐ Disassembled/Decontaminated

☐ Inventory Check Complete

☐ All Items Complete and Accounted For

Returned By: _____

Received By: _____

Date: _____

Notes: _____

WWW.BATRIK.COM TAG-IT™ TI-35L3

BROKEN / ORANGE

ITEM #	SIZE
TI-3507	3" x 5" (8 cm x 13cm)



BROKEN DO NOT USE

DATE: _____ PHONE: _____

DEPARTMENT/CONTACT: _____


ITEM: _____

CRITICAL: (必) Yes No

WWW.BATRIK.COM TAG-IT™ TI-3507

BLANK / BLACK

ITEM #	SIZE
TI-3509	3" x 5" (8 cm x 13cm)



WWW.BATRIK.COM TAG-IT™ TI-3509



QUICK TURNOVER

Tech Initials: _____

Room Number: _____

Item Number: _____

Time Needed: _____

☐ Decontam/LLD/ILD/RLD

☐ Immediate Use Sterilization

☐ Terminal Sterilization in SPD

WWW.BATRIK.COM TAG-IT™ TI-35T2

QUICK TURNOVER/ BLUE

ITEM #	SIZE
TI-35T2	3" x 5" (8 cm x 13cm)

All information tags are single-use. PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.

NEW!**BATRIK**

TAG-IT™ | SCOPE

Batrik Scope Tags are **tear-resistant, double-sided** tags that will prove indispensable throughout the reprocessing journey. Our tags make it easy to **record important information**, including processing dates, with a user-friendly feature that allows you to select the month, day, and date. Crafted from **waterproof synthetic paper**, these tags feature a **self-looping design**, making it effortless to attach them to the scope without the need for additional fastening materials. Elevate your reprocessing protocols with Batrik!

KEY FEATURES

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C) and are compatible with Gas plasma sterilization (STERRAD®100 NX Sterilization System).
- Tags include the month, day, and date selection on one side, making it convenient to mark the processing date accurately



COLOR-CODED AND CUSTOMIZABLE

- Available in seven colors
- The tags facilitate color coding for easy identification

SELF LOOPING

- The self-looping design allows for easy attachment to the scope
- Eliminating the need for additional fastening materials


SINGLE-USE AND LATEX FREE

- Batrik Scope Tags are designed for single use and are latex-free
- Promoting hygiene and reducing the risk of cross-contamination

PRODUCT SELECTION

ITEM #	COLOR	SIZE	PACKAGING
TI-ST05	Yellow	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST02	Blue	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST03	Green	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST07	Orange	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST06	Purple	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST04	Red	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG
TI-ST01	White	11.8" x 1.55" / 33.8 cm x 4 cm	100 / PKG

All information tags are single-use. PACKAGING: 1 pack of 100 tags and 1 permanent marker.



NEW!

BATRIK

TAG-IT™ | SCOPE REPROCESSING

Batrik Scope Reprocessing Tags are durable and can withstand handling during reprocessing **without tearing**, ensuring that information remains intact. These tags are vital in recording important information such as **date to reprocess on, scope model and serial number, sterilizer type, meeting sterilizer parameter conditions, alcohol flush, and hang date**. These tags have been meticulously designed to elevate **operational efficiency** and contribute to the delivery of the highest standard of patient care.




LOGGING IMPORTANT INFORMATION

- Easy for staff to indicate when the scope was processed
- Key scope, technician, and cycle-related information easily reported and identified.

SMART DESIGN

- The self-looping design allows for easy attachment to the scope
- No requirement for extra fastening materials
- Secure fastening as the loop is tightened around the instrument
- Features an integrated barcode area

TI-RST



TI-MPB



*Suitable for various information tagging purposes

PRODUCT SELECTION



ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
TI-RST	REPROCESSING SCOPE TAG	3" x 10.5" / 8 cm x 27 cm	100 / PKG
TI-MPB	MULTI-PURPOSE BLANK TAG	3" x 10.5" / 8 cm x 27 cm	100 / PKG

All information tags are single-use. PACKAGING: 1 pack of 100 tags and 1 permanent marker.


BATRIK SILICONE PIN MAT

Our Silicone Pin Mat cushions, protects, and organizes delicate surgical instruments for sterilization. The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.

- Reusable
- Autoclavable
- Latex free



FRONT SIDE
Pins keep instruments securely in place



BACK SIDE
Numerous holes and elevated nubs provide superior venting and reduce water collection

PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
13-2001	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	1 unit
13-2005	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	5 units

GOOD TO KNOW:
CUSTOM SIZES
ARE AVAILABLE
UPON REQUEST

BATRIK CORNER POCKETS™


With three models to choose from, our Corner Pockets™ portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



18-CP02



18-CP07




18-CP06

COLOUR	BLUE	ORANGE	PINK
USAGE	SINGLE-USE	REUSABLE <small>withstands repeated sterilization up to 134°C</small>	REUSABLE <small>withstands repeated sterilization up to 134°C</small>
DIMENSIONS	2.8" x 4" (7.3 cm x 10.3 cm)	2.8" x 4" (7.3 cm x 10.3 cm)	3.4" x 4.7" (8.5 cm x 12 cm)
LATEX-FREE	✓	✓	✓
NO MOISTURE RETENTION	✓	✓	✓
FABRICATION	Flexible, silicone-like plastic	Flexible, silicone-like plastic	Our most flexible plastic for easiest placement and most secure protection

SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors are made from high density foam to effectively cushion and protect the corners of surgical trays to prevent tearing. They are especially useful with heavier instrument sets.



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
18-CP02 **	Corner Pockets™ Single-Use, Blue	240 / Bag
18-EVAL2 **	Corner Pockets™ Single-Use, Blue	24 / Bag
18-CP07	Corner Pockets™ Reusable, Orange	80 / Bag
18-CORNER	Corner Pockets™ Reusable, Orange	320 / Box (4 Bags of 80 pcs)
18-EVAL7	Corner Pockets™ Reusable, Orange	24 / Bag
18-CP06	Corner Pockets™ Reusable, Pink	80 / Bag
18-CPF-02	Soft Pockets™ Single-Use Foam Corner Protectors	80 / Box (24 Boxes / Pallet)

** Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.

BATRIK

LILA'S INSTRUMENT STAIN REMOVER

Lila's special blend of naturally sourced powders provides outstanding removal of stubborn stains.

Convenient-to-use formula, you can achieve spotless results every time!

It can be used on surgical/dental instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, trays, bowls and more.

Specifically designed for use on stainless steel, chrome, nickel, brass, copper, and silver.




EASE TO USE AS 1,2,3...

1. Dampen a cloth or grab a Lila's Low Linting Dura-Towel and dab into the jar.
2. Gently rub the area that is stained.¹
3. Rinse thoroughly for a flawless shine.²

A little Goes a long Way!

Even great to tackle stubborn stains on steel pots and pans in your kitchen with ease.

Lila's - If you have used Spectra-Scrub or Integra's instrument stain remover then you will love Lila's!

**LSR-8501****17-DT1213C**

LILA'S SCRUB PAD

- Say goodbye to extra stubborn stains with Lila's Scrub Pad without leaving any metal residue.
- Simply add a small amount of Lila's Stain Remover with one of our Scrub Pad's for very large areas.
- The gentle yet effective scrub is the perfect choice for Ophthalmic Instruments.¹

LILA'S LOW-LINTING DURA-TOWEL

- Lila's Dura-Towel has long lasting strength and has a minimal lint composition.
- These single-use towels are clean room certified and are specifically designed for instrument wiping.
- Its ideal dimensions and unique attributes guarantee optimal performance, setting it apart from ordinary paper towels.

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
LSR-8501	LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid	Ø 2.8" / 70 mm	1 Jar / Bag
LSR-8512	LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid	Ø 2.8" / 70 mm	12 Jars / Box
17-DT1213C	LILA'S LOW LINTING DURA-TOWEL 1080 towels	12" x 13" / 31 cm x 33 cm Folded: 3" x 3.25" / 7.5 cm x 8.25 cm	90 / Bag 12 Bags / Box
60-46003	LILA'S SCRUB PAD	4" x 6" / 10 cm x 15 cm	3 / Bag
LSR-BUNDLE	LILA'S STAIN REMOVER + LILA'S DURA-TOWEL 3 oz (85g) plastic jar + 90 towels	Ø 2.8" / 70 mm 12" x 13" / 31 cm x 33 cm	1 Jar / Bag 90 / Bag

1 : Although Lila's is completely safe, it is not recommended for use on ophthalmic instruments. Use Lila's Scrub Pad on its own.

2: The stain remover should not be mixed with bleach and is for external use only to remove instrument stains.


BATRIK

INSTRUMENT PROTECTORS


Our Instrument Protectors protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization*
- Multiple color-coded sizes
- Latex Free
- Non-Sterile

REGULAR / CLEAR



PERFORATED / CLEAR



REGULAR / TINTED



PERFORATED / TINTED



REGULAR / SOLID COLOR



PERFORATED / SOLID COLOR



*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy in a healthcare setting

PRODUCT SELECTION

	ITEM # REGULAR	ITEM # PERFORATED	COLOR	DIMENSIONS
 RC-07  PC-07	REGULAR & PERFORATED / CLEAR			
	RC-01	PC-01	Clear	0.062" x 0.75" (1.57 mm x 19 mm)
	RC-02	PC-02	Clear	0.075" x 0.75" (2 mm x 19 mm)
	RC-03	PC-03	Clear	0.109" x 0.75" (2.77 mm x 19 mm)
	RC-04	PC-04	Clear	0.125" x 1" (3.18 mm x 25.4 mm)
	RC-05	PC-05	Clear	0.187" x 1" (4.75 mm x 25.4 mm)
	RC-06	PC-06	Clear	0.406" x 0.75" (10.3 mm x 19 mm)
	RC-07	PC-07	Clear	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RC-08	PC-08	Clear	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RC-09	PC-09	Clear	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RC-10	PC-10	Clear	Assorted
PACKAGING: Bag of 100				
 RT-07  PT-07	REGULAR & PERFORATED / TINTED			
	RT-01	PT-01	Tinted White	0.062" x 0.75" (1.57 mm x 19 mm)
	RT-02	PT-02	Tinted Blue	0.075" x 0.75" (2 mm x 19 mm)
	RT-03	PT-03	Tinted Green	0.109" x 0.75" (2.77 mm x 19 mm)
	RT-04	PT-04	Tinted Red	0.125" x 1" (3.18 mm x 25.4 mm)
	RT-05	PT-05	Tinted Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
	RT-06	PT-06	Tinted Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
	RT-07	PT-07	Tinted Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RT-08	PT-08	Tinted Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RT-09	PT-09	Tinted Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RT-10	PT-10	Tinted Assorted	Assorted
PACKAGING: Bag of 100				
 RS-07  PS-07	REGULAR & PERFORATED / SOLID			
	RS-01	PS-01	Solid White	0.062" x 0.75" (1.57 mm x 19 mm)
	RS-02	PS-02	Solid Blue	0.075" x 0.75" (2 mm x 19 mm)
	RS-03	PS-03	Solid Green	0.109" x 0.75" (2.77 mm x 19 mm)
	RS-04	PS-04	Solid Red	0.125" x 1" (3.18 mm x 25.4 mm)
	RS-05	PS-05	Solid Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
	RS-06	PS-06	Solid Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
	RS-07	PS-07	Solid Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RS-08	PS-08	Solid Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RS-09	PS-09	Solid Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RS-10	PS-10	Solid Assorted	Assorted
PACKAGING: Bag of 100				

ITEM #	DIMENSIONS
OP-11	0.25" x 0.75" (6.35 mm x 19.05 mm)
OP-12	0.375" x 1" (9.53 mm x 25.4 mm)
OP-13	0.5" x 1" (12.7 mm x 25.4 mm)
OP-14	0.75" x 1" (19.05 mm x 25.4 mm)
OP-15	1" x 1" (25.4 mm x 25.4 mm)
OP-16	1.25" x 1.25" (31.75 mm x 31.75 mm)
OP-17	1.5" x 1.5" (38.1 mm x 38.1 mm)
OP-18	1.75" x 1.25" (44.45 mm x 31.75 mm)
OP-20	Assorted

PACKAGING: Bag of 50

ITEM #	DIMENSIONS	PACKAGING
RGS-01	1.5" (38 mm)	55 / Bag
RGS-02	2" (51 mm)	45 / Bag
RGS-03	3" (76.2 mm)	40 / Bag


ITEM #	COLOR	SIZE
RGM-01	Pink	Small
RGM-02	Clear	Small
RGM-03	Green	Medium
RGM-04	Clear	Medium

PACKAGING: Bag of 50

ITEM #	COLOR	DIMENSIONS
DG-01	Clear	0.0625" x 0.75" (1.59 mm x 19.05 mm)
DG-02	Blue	0.078" x 0.75" (1.98 mm x 19.05 mm)
DG-03	Green	0.109" x 0.75" (2.78 mm x 19.05 mm)
DG-04	Red	0.125" x 1" (3.18 mm x 25.4 mm)
DG-05	Yellow	0.187" x 1" (4.76 mm x 25.4 mm)
DG-07	Orange	0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm)
DG-10	Assorted	Assorted

PACKAGING: Bag of 50


OSTEOTOME PROTECTORS




SINGLE-TIP RETRACTOR PROTECTORS



WEITLANER RETRACTOR PROTECTORS



PERFORATED DOUBLE PROTECTORS



BATRIK PrAmLi Power³ LUBE™ Instrument Lubricant

Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving. PrAmLi Power³ LUBE™ is uniquely formulated with an all natural corrosion inhibitor.

- Forms a protective coating which can extend the instrument's, useful life
- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion

TIP:
FOR MILK BATH SOLUTION,
DILUTE 1 PART SOLUTION
IN 6 PARTS WATER

Formulated with
PrAmLi Power³ Enzymes™



PRODUCT SELECTION


ITEM #	DILUTION LUBRICANT / WATER	PACKAGING
SS-21-6001	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	1 x 3.8 L (1 Gallon)
SS-21-6004	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	4 x 3.8 L (4 Gallons)
SS-21-6025	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)
SS-21-6032	No dilution required	1 x 946 ml (32 oz) Spray bottle
SS-21-6012	No dilution required	12 X 946 ml (32 oz) Spray Bottles

BATRIK S.T.I.C.S.

Our Simple Tracking Identification Color System, or S.T.I.C.S., is an easy to use, latex free, tool that allows you to quickly identify your instruments by department and/or procedure. The tape adheres securely and every tape sheet comes with a protective cover to ensure it lasts as long as you need.

Each 8.5" x 11" (216 mm x 279 mm) sheet contains the following adhesive strips:

- 1" (25.4 mm) x 220
- 1.5" (38.1 mm) x 44
- 2" (50.8 mm) x 44
- Total tape per sheet 374" (9500 mm)



PRODUCT SELECTION

SOLID COLORS	
ITEM #	COLOR
T-850-2-01	White
T-850-2-02	Royal Blue
T-850-2-03	Kelly Green
T-850-2-04	Red
T-850-2-05	Yellow
T-850-2-06	Purple
T-850-2-07	Orange
T-850-2-08	Brown
T-850-2-09	Black
T-850-2-10	Gray
T-850-2-11	Lavender
T-850-2-12	Beige
T-850-2-13	Light Blue
T-850-2-14	Lime Green
T-850-2-15	Pink
T-850-2-17	Plum
T-850-2-18	Mint Green

BLACK STRIPES	
ITEM #	COLOR
T-850-3-01	Black/White
T-850-3-03	Black/Kelly Green
T-850-3-04	Black/Red
T-850-3-05	Black/Yellow
T-850-3-07	Black/Orange
T-850-3-10	Black/Gray
T-850-3-11	Black/Lavender
T-850-3-12	Black/Beige
T-850-3-13	Black/Light Blue
T-850-3-14	Black/Lime Green
T-850-3-15	Black/Pink
T-850-3-17	Black/Plum

TRI-COLOR	
ITEM #	COLOR
T-850-5-51	Red/White/Blue
T-850-5-54	Red/Black/Green

SPECIALTY	
ITEM #	COLOR
T-850-7-04	Yellow/Red


WHITE STRIPES	
ITEM #	COLOR
T-850-4-02	White/Royal Blue
T-850-4-03	White/Kelly Green
T-850-4-04	White/Red
T-850-4-05	White/Yellow
T-850-4-07	White/Orange
T-850-4-08	White/Brown
T-850-4-11	White/Lavender
T-850-4-13	White/Light Blue
T-850-4-14	White/Lime Green
T-850-4-15	White/Pink
T-850-4-17	White/Plum

FOUR-COLOR	
ITEM #	COLOR
T-850-2-31	Red/Yellow/ Purple/Orange
T-850-2-32	Blue/White/ Black/Green

BATRIK ID TAPE



AVAILABLE WIDTHS



A simple solution for marking surgical instruments with the aid of color-coding, our ID Tape is durable and will not chip, flake, or peel. It can be used in autoclaves, withstanding repeated exposure to sterilization. Independent validation tests demonstrated the efficiency of sterilization of marked instruments, even underneath the BATRIK ID Tape.



Quickly identifies instruments by department or procedure



Easy to apply



Withstands repeated exposure to sterilization



Has more per tape roll (300" / 7620 mm)



Available in 9 colors and 2 widths



Conveniently packaged in a clear box



Non-Sterile

PRODUCT SELECTION

ITEM # 0.125" / 318 MM	ITEM# 0.25" / 635 MM	COLOR
T-125-01	T-250-01	White
T-125-02	T-250-02	Blue
T-125-03	T-250-03	Green
T-125-04	T-250-04	Red
T-125-05	T-250-05	Yellow
T-125-06	T-250-06	Purple
T-125-07	T-250-07	Orange
T-125-08	T-250-08	Brown
T-125-09	T-250-09	Black

BATRIK

TAPE-B-GONE™ TAPE REMOVER

Tape-B-Gone™ is a uniquely formulated and ready-to-use proprietary formula that easily removes adhesives from medical equipment and instruments. Unlike many other brands, Tape-B-Gone™ does not contain highly dangerous or caustic chemicals, is not flammable or combustible, and is safe for use on metal, glass, plastic, and rubber without any special equipment or ventilation.



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
T-08RMV	8 oz (237 ml) bottle	1 Unit
T-32RMV	8 oz (237 ml) bottle	4 / Pack
T-SAMP-RMV	0.5 oz (15 ml) bottle	1 Unit

BATRIK


PADLOCKS SECURITY TAGS

Our Padlock Security Tags for instrument containers and general hospital applications are specifically designed for the healthcare market. Our material selection, fabrication process, and engineering ensure consistent quality and high temperature compliance with the various steam sterilization temperatures in worldwide hospital applications.


Available with optional reactive ink indicator dots that provide visual sterility indication. Light blue dots turn brown/dark brown when steam sterilized.

AVAILABLE IN THREE FORMATS


700-SERIES



800-SERIES



900-SERIES



Made of polypropylene



Latex and DEHP Free



Will not melt or soften with any established steam sterilization protocol



Fully open design (45°) facilitates insertion into instrument containers



Easy twist breaking design



Large base for easier handling



Available in different size



Custom Colors and OEM available upon request

PRINTED & CUSTOM NUMBERED PADLOCKS

Available in 700, 800 & 900 SERIES

PRODUCT SELECTION

COLOR	PRINTED WORD
PRINTED PADLOCKS	
WHITE	CLEAN
BLUE	STERILE
RED	REPAIR
YELLOW	SHARPEN
YELLOW	LOANER
BLACK	DISCARD
CUSTOM NUMBERED	
ALL COLORS	up to 8 numbers



CUSTOM NUMBERED
Custom printed padlocks (up to 8 characters)



PRINTED
6 different printed labels available

NOTE: Printed & Custom Numbered Padlocks are available in 700, 800, and 900 Series and can be ordered in 100, 200, and 500 unit quantities. Please contact us for ordering details.

INDICATOR DOTS



REGULAR



700-SERIES PADLOCKS

The 700 series padlocks are perfect for "Aesculap" containers and those with a short, large opening.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701	15-711	15-710
BLUE	15-702	15-712	15-720
GREEN	15-703	15-713	15-730
RED	15-704	15-714	15-740
YELLOW	15-705	15-715	15-750
PINK	15-706	15-716	15-760
ORANGE	15-707	15-717	15-770
BLACK	15-709	15-719	15-790
LIGHT BLUE	15-722	15-732	15-742

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701-IND	15-711-IND	15-710-IND
BLUE	15-702-IND	15-712-IND	15-720-IND
GREEN	15-703-IND	15-713-IND	15-730-IND
RED	15-704-IND	15-714-IND	15-740-IND
YELLOW	15-705-IND	15-715-IND	15-750-IND
PINK	15-706-IND	15-716-IND	15-760-IND
ORANGE	15-707-IND	15-717-IND	15-770-IND
BLACK	15-709-IND	15-719-IND	15-790-IND
LIGHT BLUE	15-722-IND	15-732-IND	15-742-IND

Double locking hooks can be inserted in a large guide hole



The **grooves** on the padlock facilitate handling



INDICATOR
DOTS

REGULAR



800-SERIES PADLOCKS

The 800 series padlocks have a longer arm (32 mm) and have a single, thinner hook than the 700 series. They are perfect for almost any type of container, even those with long and thin channels.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801	15-811	15-810
BLUE	15-802	15-812	15-820
GREEN	15-803	15-813	15-830
RED	15-804	15-814	15-840
YELLOW	15-805	15-815	15-850
PINK	15-806	15-816	15-860
ORANGE	15-807	15-817	15-870
BLACK	15-809	15-819	15-890
LIGHT BLUE	15-822	15-832	15-842

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801-IND	15-811-IND	15-810-IND
BLUE	15-802-IND	15-812-IND	15-820-IND
GREEN	15-803-IND	15-813-IND	15-830-IND
RED	15-804-IND	15-814-IND	15-840-IND
YELLOW	15-805-IND	15-815-IND	15-850-IND
PINK	15-806-IND	15-816-IND	15-860-IND
ORANGE	15-807-IND	15-817-IND	15-870-IND
BLACK	15-809-IND	15-819-IND	15-890-IND
LIGHT BLUE	15-822-IND	15-832-IND	15-842-IND

Single locking hook
can be inserted
in a large guide hole



The locking system
allows for **easy closure**



INDICATOR DOTS



REGULAR



900-SERIES PADLOCKS

The 900 series padlocks are identical to the 800 series, but with a larger rectangular base which allows for the placement of a 20 mm x 30 mm tracking label.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901	15-911	15-910
BLUE	15-902	15-912	15-920
GREEN	15-903	15-913	15-930
RED	15-904	15-914	15-940
YELLOW	15-905	15-915	15-950
PINK	15-906	15-916	15-960
ORANGE	15-907	15-917	15-970
BLACK	15-909	15-919	15-990
LIGHT BLUE	15-922	15-932	15-942

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901-IND	15-911-IND	15-910-IND
BLUE	15-902-IND	15-912-IND	15-920-IND
GREEN	15-903-IND	15-913-IND	15-930-IND
RED	15-904-IND	15-914-IND	15-940-IND
YELLOW	15-905-IND	15-915-IND	15-950-IND
PINK	15-906-IND	15-916-IND	15-960-IND
ORANGE	15-907-IND	15-917-IND	15-970-IND
BLACK	15-909-IND	15-919-IND	15-990-IND
LIGHT BLUE	15-922-IND	15-932-IND	15-942-IND

Single locking hook
for use with smaller
diameters



Ideal for adding labels/
indicators indicators or
markings to the large
rectangular base



BATRIK SAVEFACE™ FACE SHIELDS

Our Batrik SaveFace™ Face Shields (Visors) are lightweight and allow for good peripheral vision while acting as a physical barrier against splashes. Can be used for extended periods due to comfort-focused design and superior quality materials.

Choose amongst non-draped and draped models:

- Our classic **non-draped** models close the gap between the shield and the chest bone when head is tilted forward for added protection.
- Our **draped** models offer full-length coverage extending from the top of the eyebrow to below the chest, covering the anterior neck area as per CDC recommendations.

SF-100

SF-240



SF-180

SF-8240



HEADPIECE

- Self-shaping foam headpieces mold and cushion forehead for a comfortable and secure fit (per CDC recommendations)
- Shielded air vents maximize airflow for enhanced comfort and prevent fogging

VISOR

- Excellent Distortion-Free Optical Clarity with flat viewing area reduces distortion and glare
- Features anti-fog treatment for clear view
- Wide visors provide peripheral protection to help prevent splashes from going into the eyes
- High Quality PET
 - SF-100/SF-180: 7.5 mil (0.19 mm) thick
 - SF-240/SF-8240: 9.84 mil (0.25 mm) thick

DRAPE

- Our patent-pending drape models come equipped with integrated ANSI/AAMI PB70 Level 3 compliant splash drape barriers for added protection



INCLUDES CLOSING SYSTEM TO
REDUCE RISK OF FLUID CONTACT

PATENT PENDING

ASSEMBLED IN CANADA 

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-100	SAVEFACE™ no drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm	9" (23 cm)	100 / Box
SF-180	SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm Drape Length: 6" / 15 cm	15" (38 cm)	100 / Box
SF-240	BATRIK SAVEFACE™ no drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm	9" (23 cm)	240 / Case (8 Boxes of 30 Pcs)
SF-8240	BATRIK SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm Drape Length: 6" / 15 cm	15" (38 cm)	240 / Case (8 Boxes of 30 Pcs)

BATRIK SHIELD™ FACE SHIELD

A cost-effective alternative to our premium SAVEFACE™ models, the BATRIK Shield™ provides excellent splash protection. An ideal safety addition to all pandemic kits.

- Self-shaping foam head piece molds comfortably to forehead
- Shields can be used for longer periods due to comfort and improved materials
- Optical shield with anti-fog and anti-glare provide a clearer viewing area
- Soft plastic fasteners adhere the elastic band to the shield for additional comfort and easy use
- Latex free



PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-096	BATRIK Shield™ no drape	Shield Length: 9" / 23 cm*	9" (23 cm)*	96 / Box
SF-896	BATRIK Shield™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm* Drape Length: 7" / 17.8 cm	15.9" / 40.3 cm*	96 / Box

*Shield product tolerances: 8.9" - 9.5" / 22.5 cm - 23.2 cm

BATRIK SCRUB PAD

Our Scrub Pad is ideal for pre-cleaning your surgical instruments without leaving any metal residue that often occurs when using steel wool or metal sponges. The abrasiveness is perfectly balanced to provide you with efficient cleaning without damaging your instruments. Its fine fibers allow the user to polish, clean, and debur all types of material surfaces.

- Abrasive Material: Aluminum Oxide
- Grit: 400
- Grade: Very Fine
- Style: Flexible
- Pad material: Non-Woven
- Color: Maroon
- Type: Hand Pad



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
60-46050	Scrub Pad	4" x 6" / 10 cm x 15 cm	150 / Box (50 bags of 3 units each)
60-46003	Scrub Pad	4" x 6" / 10 cm x 15 cm	3 / Bag

Your First ™

ENZYMATIC & NON-ENZYMATIC CLEANERS

Our Enzymatic and Non-Enzymatic Cleaners are specially formulated to remove organic material from medical equipment and instrumentation.

Our **Enzymatic Cleaner** formats are formulated with PrAmLi Power3 Enzymes™, our proprietary unique blend of enzymes, and a surfactant to efficiently remove all proteins, blood, fatty tissue, and more. They also remove stains from surgical scrubs, linens, and medical clothing.

Our **Non-Enzymatic Cleaner** formats are formulated with a free-rinsing, neutral pH detergent blend for effective cleaning.



Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning*

PRODUCT SELECTION

ENZYMATIC CLEANER	NON-ENZYMATIC CLEANER	DILUTION DETERGENT / WATER	PACKAGING
LOW FOAM			
SS-21-1001	SS-721-1001	4 to 8 ml / L (0.5 to 1 oz / Gal)	1 x 3.8 L (1 Gal)
SS-21-1004	SS-721-1004	4 to 8 ml / L (0.5 to 1 oz / Gal)	4 x 3.8 L (4 Gal)
MILD FOAM			
SS-21-2001	SS-721-2001	4 to 8 ml / L (0.5 to 1 oz / Gal)	1 x 3.8 L (1 Gal)
SS-21-2004	SS-721-2004	4 to 8 ml / L (0.5 to 1 oz / Gal)	4 x 3.8 L (4 Gal)
SUPER CONCENTRATE			
SS-21-3001	SS-721-3001	1 to 2 ml / L (0.125 to 0.25 oz / Gal)	1 x 3.8 L (1 Gal)
SS-21-3004	SS-721-3004	1 to 2 ml / L (0.125 to 0.25 oz / Gal)	4 x 3.8 L (4 Gal)

NOTE: Our Low Foam & Super Concentrate Detergent formulas are available in a 2.5 Gallon format based on a minimum order quantity. Contact us for more information.

GOOD TO KNOW:

OUR LOW FOAM & SUPER CONCENTRATE DETERGENTS (ENZYMATIC & NON-ENZYMATIC) MAY BE USED IN WASHERS / ULTRASONIC CLEANERS.


WE RECOMMEND CONFIRMING SUGGESTED DILUTION RATIO AND WATER TEMPERATURE WITH THE MANUFACTURER OF YOUR WASHER / ULTRASONIC CLEANER AS MACHINES / SETTINGS VARY.

*Caution: We recommend consulting your hospital pre-cleaning/cleaning protocols for various equipment/instrument cleaning detergent preferences.

Your First ™

ENVIRO-SACHET™

ENZYMATIC GEL & NON-ENZYMATIC CLEANER

Your First ™ Enviro-Sachet™, our unique cleaning solutions now available in Enzymatic & Non-Enzymatic formats.



- Our Enzymatic Gel format is formulated with the same triple action **PrAmLi Power³ Enzymes™** that you know and trust.
- Our Non-Enzymatic Cleaner format is formulated with a free-rinsing, neutral pH detergent blend.

Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning.

PRODUCT SELECTION

ENZYMATIC GEL	NON-ENZYMATIC CLEANER	DESCRIPTION	PACKAGING
17-166	717-166	1 Dry Tubular Sponge & Enviro-Sachet™ set	25 / Box 4 Boxes / Case
17-266	717-266	1 Dry Flat Sponge & Enviro-Sachet™ set	50 / Box 4 Boxes / Case
17-177	717-177	1 Dry Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit
17-277	717-277	1 Dry Flat Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit
17-EZSAMP	717-EZSAMP	SAMPLE KIT INCLUDES: 3 Dry Flat Sponges with Enviro-Sachet™ sets 1 Dry Tubular Sponge with Enviro-Sachet™ set 1 Reusable Stainless Steel Bowl	1 Kit
H-EZ500		Reusable Stainless Steel Bowl 4.3" x 5.3" x 3" / 11 cm x 13.5 cm x 7.5 cm	1 Unit

Your First ™

FOAMING SPRAY

READY-TO-USE ENZYME FOAMING SPRAY

Our Foaming Spray is ideal for instrument pre-cleaning as the thick foam spray enzymatically breaks down all bioburden to allow for effective cleaning and sterilization.

Formulated with **PrAmLi Power³ Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



INSTRUCTIONS



PRODUCT SELECTION

ITEM #	PACKAGING
SS-21-4001	1 x 946 ml (32 oz) Spray Bottle
SS-21-4012	12 x 946 ml (32 oz) Spray Bottles

Important: In addition to our enzymatic formulation, our spray bottle contains as much as 25% more solution when compared to the competition.

Your First ™

F1RST WIPES

Our cleanroom wipes are designed to do the first wipe-down of an endoscope after a procedure. They are made from hydroentangled non-woven material composed of 45% polyester and 55% cellulose. These wipes are manufactured and packaged in a cleanroom environment.



- Can be used for a cleanroom environment of ISO Class 5 through 8
- Production process guarantees extremely low particle generation and chemical extractables
- Compatible with virtually any solvent or solution for cleaning and wiping (such as **PrAmLi Power³ Enzymes™**)
- Double-bag packaging
- Available in two versatile sizes


PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
17-A0909	Non-woven Cleanroom Wipes	9" x 9" (23 cm x 23 cm)	300 / Bag 14 Bags / Case
17-A1212	Non-woven Cleanroom Wipes	12" x 12" (31 cm x 31 cm)	150 / Bag 10 Bags / Case

Your First ™

DRY SPONGES

Our sponges are specifically fabricated for use with all medical instrumentation. Unlike most competitor sponges, our sponges are white to better evidence the effectiveness of the cleaning process. Dark colored sponges easily hide contaminants, but with our white sponges you can visually evaluate the pre-clean efficacy.



OUR WHITE SPONGES MAY APPEAR SLIGHTLY YELLOWED. THIS IS ONLY DUE TO UV EXPOSURE, IS NOT A REFLECTION OF ITS AGE.

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
17-100	DRY TUBULAR SPONGE White Dispenser Box	1.06" ø x 4.25" (2.7 cm ø x 10.8 cm)	25 / Box (4 Boxes / Case)
17-101	DRY TUBULAR SPONGE Bulk	1.06" ø x 4.25" (2.7 cm ø x 10.8 cm)	560 / Case
17-200	DRY FLAT SPONGE White Dispenser Box	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	50 / Box (4 Boxes / Case)
17-201	DRY FLAT SPONGE Bulk	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	1,000 / Case

Your First ™

ENZYMATIC SPONGES & BEDSIDE SPONGE KITS



DISPENSER BOX

Our Enzymatic Bedside Sponge Kits feature sponges all pre-treated with **PrAmLi Power³ Enzymes™** to provide the most effective pre-clean or surface clean solution. This unique three-enzyme blend provides you with confident and effective cleaning. The tubular sponge is ideal for all endoscopes.



FREESTANDING POUCH



FLAT SPONGE



TUBULAR SPONGE

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
17-111	PRE-TREATED TUBULAR SPONGE Add 100-500 ml of water to activate detergent solution.	25 / Box (4 Boxes / Case)
17-211	PRE-TREATED FLAT SPONGE Add 100-500 ml of water to activate detergent solution.	50 / Box (4 Boxes / Case)
17-122	PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add up to 500 ml of water.	15 / Box (4 Boxes / Case)
17-222	FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add up to 500 ml of water.	25 / Box (4 Boxes / Case)
17-133	READY TO USE TUBULAR SPONGE Pre-treated in 500 ml of an enzymatic solution in freestanding pouch.	24 / Case
17-233	READY TO USE FLAT SPONGE Comes bathing in 500 ml of enzymatic solution.	24 / Case
17-144	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes in a freestanding pouch soaked in 100 ml of enzymatic solution.	12 / Box (4 Boxes / Case)
17-244	READY TO USE DUAL PURPOSE FLAT SPONGE Comes in a freestanding pouch soaked in 100 ml of enzymatic solution.	12 / Box (4 Boxes / Case)
17-155	READY TO USE TUBULAR SPONGE Comes bathing in 220 ml of enzymatic solution.	24 / Case
17-255	READY TO USE FLAT SPONGE Comes bathing in 220 ml of enzymatic solution.	24 / Case

NON-ENZYMATIC DETERGENT BEDSIDE SPONGE KITS



DISPENSER BOX




FREESTANDING POUCH



FLAT SPONGE



TUBULAR SPONGE

Our Non-Enzymatic Bedside Sponge Kits are all pre-treated with a free-rinsing, neutral pH detergent blend, for an effective instrument pre-cleaning and surface clean solution. This kit is **Your First ™** in the medical instrument cleaning process.

PRODUCT SELECTION


ITEM #	DESCRIPTION	PACKAGING
717-122	PRE-TREATED TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	15 / Box (4 Boxes / Case)
717-222	PRE-TREATED FLAT SPONGE Comes in a freestanding pouch to use as a basin. dd 500 ml of water.	25 / Box (4 Boxes / Case)
717-133	READY TO USE TUBULAR SPONGE Comes bathing in 500 ml of detergent solution.	24 / Case
717-233	READY TO USE FLAT SPONGE Comes bathing in 500 ml of detergent solution.	24 / Case
717-144	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes in a freestanding pouch soaked in 100 ml of detergent solution.	12 / Box (4 Boxes / Case)
717-244	READY TO USE DUAL PURPOSE FLAT SPONGE Comes in a freestanding pouch soaked in 100 ml of detergent solution.	12 / Box (4 Boxes / Case)
717-155	READY TO USE TUBULAR SPONGE Comes bathing in 220 ml of detergent solution.	24 / Case
717-255	READY TO USE FLAT SPONGE Comes bathing in 220 ml of detergent solution.	24 / Case

Your First ™

TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

Thoughtfully designed to eliminate touch-points during endoscope transport from point-of-care to the reprocessing area, and back, focused on safety of patients and medical staff.

- Made of recyclable LDPE material known to be quite flexible, highly resilient, and tough
- Latex and DEHP-free
- Excellent resistance to dilute and concentrated acids, alcohols, bases and esters*
- Good resistance to aldehydes, ketones, and vegetable-based oils*
- Can withstand high temperatures of up to 90 °C (194°F)
- Tamper-evident seal for added safety, also aids in keeping the instrument moist to avoid drying of bio-debris
- Clearly marked bags with color-coded visuals facilitates immediate identification of clean vs. used instruments
- Users can write important information directly on the bag in spaces marked with visual icons
- Available in two convenient sizes


GOOD TO KNOW: USING **Your First ™** SEALABLE TRANSPORT BAGS REDUCES USER ERROR, PREVENTS CROSS-CONTAMINATION, AND PROMOTES PATIENT AND MEDICAL STAFF SAFETY.

PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
17-HAP24	Clean Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case
17-SAD24	Soiled Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case
17-HAP30	Clean Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case
17-SAD30	Soiled Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case

*Reference on file

CLEAN INSTRUMENTS BAG



SOILED INSTRUMENTS BAG



MADE IN CANADA 

BATRIK

ARC™ WIRE STYLE BRUSHES

SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

BRUSH SELECTION TIPS:

INSTRUMENT SPECIFICATIONS

What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

What is the inner diameter of the instrument?

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that are too large will bend too much and not clean effectively.

BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.


In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At BATRIK, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

CUSTOMIZABLE WIRE BRUSHES


A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000mm for Wire Brushes on plastic semi-rigid sheaths, rigid brushes from 5 to 20mm in diameter.
- No limit in length for Wire Brushes on flexible plastic sheaths, flexible brushes from 2-5mm in diameter

RIGID



FLEXIBLE



12" / 30 CM SINGLE-END WIRE BRUSH - DELRIN PLASTIC

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-31502-S	41-31502	White	12.401" / 31.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-32503-S	41-32503	Yellow	12.795" / 32.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-32504-S	411-32504	Violet	12.795" / 32.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-32505-S	411-32505	Grey	12.795" / 32.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-35005-S	41-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35006-S	41-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35008-S	41-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35010-S	41-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35012-S	41-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35015-S	41-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-35020-S	41-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING Single-Use: Bag of 20 / Reusable: Bag of 3

18" / 45 CM SINGLE-END WIRE BRUSH - DELRIN PLASTIC

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-46502-S	41-46502	White	18.307" / 46.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-47503-S	41-47503	Yellow	18.7" / 47.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	41-47505-S	N/A	White	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
	411-47504-S	411-47504	Violet	18.7" / 47.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-47505-S	411-47505	Grey	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-50005-S	41-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50006-S	41-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50008-S	41-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50010-S	41-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50012-S	41-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50015-S	41-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-50020-S	41-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

24" / 60 CM SINGLE-END WIRE BRUSH - DELRIN PLASTIC

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-61502-S	41-61502	White	24.212" / 61.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1mm
	41-62503-S	41-62503	Yellow	24.6" / 62.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1mm
	41-62505-S	N/A	White	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75mm
	411-62504-S	411-62504	Violet	24.6" / 62.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55mm
	411-62505-S	411-62505	Grey	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75mm
RIGID BRUSH & SHEATH	41-65005-S	41-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3mm
	41-65006-S	41-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3mm
	41-65008-S	41-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65010-S	41-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65012-S	41-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65015-S	41-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-65020-S	41-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3							

30" / 75 CM SINGLE-END WIRE BRUSH - DELRIN PLASTIC

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-76502-S	41-76502	White	30.118" / 76.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-77503-S	41-77503	Yellow	30.511" / 77.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-77504-S	411-77504	Violet	30.511" / 77.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-77505-S	411-77505	Grey	30.511" / 77.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-80005-S	41-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80006-S	41-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80008-S	41-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80010-S	41-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80012-S	41-80012	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80015-S	41-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-80020-S	41-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3							

FLEXIBLE 2MM




FLEXIBLE 3MM



FLEXIBLE 4MM



FLEXIBLE 5MM



FLEXIBLE

Flexible brush diameter (Ø) ranging from 2 mm to 5 mm

RIGID 20MM



RIGID 15MM



RIGID 12MM



RIGID 10MM



RIGID 8MM



RIGID 6MM



RIGID 5MM



RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

BRUSH SELECTION SYSTEM (Actual Size)

2mm
0.078"



3 mm
0.118"



4 mm
0.157"



5 mm
0.196"



6 mm
0.236"



8 mm
0.314"



10 mm
0.394"



12 mm
0.472"



15 mm
0.591"



20 mm
0.787"



BATRIK

REUSABLE DOUBLE ARC™ WIRE BRUSHES

ASYMMETRIC

14" / 35 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35056	White	13.779" / 35 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35112	Yellow	13.779" / 35 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35152	Green	13.779" / 35 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35810	Blue	13.779" / 35 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm

PACKAGING: Bag of 3

20" / 50 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50056	White	19.685" / 50 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50112	Yellow	19.685" / 50 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50152	Green	19.685" / 50 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50810	Blue	19.685" / 50 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm

PACKAGING: Bag of 3

26" / 65 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65056	White	25.59" / 65 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65112	Yellow	25.59" / 65 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65152	Green	25.59" / 65 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65810	Blue	25.59" / 65 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm

PACKAGING: Bag of 3

31" / 80 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80056	White	31.496" / 80 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80112	Yellow	31.496" / 80 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80152	Green	31.496" / 80 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80810	Blue	31.496" / 80 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm

PACKAGING: Bag of 3

BATRIK

REUSABLE DOUBLE ARC™

WIRE BRUSHES

SYMMETRIC

14" / 35 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-80012	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

ASYMMETRIC BRUSHES

5 MM

6 MM



8 MM

10 MM



10 MM

12.5 MM



15 MM

20 MM




BATRIK

REUSABLE IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of Nylon Twisted Wire Brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED
WIRE BRUSHES FEATURE
A PROTECTIVE BLUE TIP
FOR ADDED SAFETY.



10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

ITEM #	SIZE	BRUSH LENGTH	BRUSH (Ø)
21-10039	3 FR	0.5" / 12.7 mm	0.039" / 1 mm
21-10067	5 FR	0.5" / 12.7 mm	0.065" / 1.67 mm
20-10059	5 FR	1" / 25.4 mm	0.059" / 1.50 mm
20-10079	6 FR	1" / 25.4 mm	0.079" / 2 mm
20-10090	7 FR	1" / 25.4 mm	0.090" / 2.3 mm
20-10106	8 FR	1" / 25.4 mm	0.106" / 2.7 mm
20-10118	9 FR	1" / 25.4 mm	0.118" / 3 mm
20-10129	10 FR	1.5" / 38 mm	0.129" / 3.3 mm
20-10146	11 FR	1.5" / 38 mm	0.0146" / 3.7 mm
20-10157	12 FR	1.5" / 38 mm	0.157" / 4 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

12" / 305 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
22-12078	2" / 50.8 mm	0.078" / 2 mm
22-12118	2" / 50.8 mm	0.118" / 3 mm
22-12157	3" / 76 mm	0.157" / 4 mm
22-12197	3" / 76 mm	0.197" / 5 mm
22-12236	3" / 76 mm	0.236" / 6 mm
22-12393	3" / 76 mm	0.394" / 10 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "- 50".

16" / 406 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
23-16197	2.5" / 63.5 mm	0.197" / 5 mm
23-16236	2.5" / 63.5 mm	0.236" / 6 mm
23-16276	3" / 76 mm	0.276" / 7 mm
23-16394	3" / 76 mm	0.394" / 10 mm
23-16472	3" / 76 mm	0.472" / 12 mm
23-16591	3" / 76 mm	0.591" / 15 mm
23-16787	3" / 76 mm	0.787" / 20 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

18" / 457 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
19-18010	0.47" / 12 mm	0.039" / 1 mm
19-18015	0.5" / 12 mm	0.06" / 1.6 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

24" / 610 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
24-24098	2.5" / 63.5 mm	0.098" / 2.5 mm
24-24197	2.5" / 63.5 mm	0.197" / 5 mm
24-24276	3" / 76 mm	0.276" / 7 mm
24-24394	3" / 76 mm	0.394" / 10 mm
24-24472	3" / 76 mm	0.472" / 12 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".



BRONCHO BRUSHES

ITEM #	DESCRIPTION	BRUSH LENGTH	BRUSH (Ø) INCHES / MM
26-06175	Small	6" / 152 mm	1.75" / 44.5 mm
26-08325	Large	8" / 203 mm	3.25" / 82.6 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

TRACHEOSTOMY BRUSHES

ITEM #	DESCRIPTION	TOTAL LENGTH	BRUSH LENGTH	BRUSH (Ø)
25-08500	Large 6-12.3 3 / 4"	8" / 203 mm	3.75" / 95 mm	0.5" / 12.7 mm
25-37550	Small 000-5.1 7 / 8"	4" / 93 mm	1.88" / 48 mm	0.5" / 12.7 mm

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".

BATRIK STAINLESS STEEL ARC™ MICRO BRUSHES

Our Stainless Steel ARC™ Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable*
- Available in two lengths

TIP: USING THE PULL-THROUGH
TECHNIQUE IS HIGHLY
RECOMMENDED WITH THESE
BRUSHES. SIMPLY INSERT THE
BLUNT, NON-BRUSH END FIRST.



PRODUCT SELECTION

ITEM #	TOTAL LENGTH	BRUSH SIZE (Ø)	BRUSH LENGTH	SHEATH (Ø)	PACKAGING
35-03015	12" / 30 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube
35-04515	18" / 45 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube

*Although all materials used in the fabrication of these brushes can be re-processed repeatedly, due to the delicate nature of these brushes re-use is conditional on their handling and application.

BATRIK

INSTRUMENT CLEANING FLAT BRUSHES

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40919



50-40921



50-40920



50-40918



50-40922



50-40917

PRODUCT SELECTION

DELTRIN PLASTIC

ITEM #	BRUSH TYPE	DESCRIPTION	COLOR	BRUSH LENGTH	BRUSH SIZE	PACKAGING
50-40917	Extra Rigid	Reusable Autoclavable Flat Brush	Blue	220 mm	75 mm x 20 mm x 10 mm	Bag of 3
50-40918	Rigid	Reusable Autoclavable Flat Brush	Green	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40919	Flexible	Reusable Autoclavable Flat Brush	Yellow	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40920	Extra Flexible	Reusable Autoclavable Flat Brush	Red	220 mm	75 mm x 20 mm x 25 mm	Bag of 3
50-40921	Flexible	Single-Use Flat Brush	Orange	220 mm	60 mm x 20 mm x 15 mm	Bag of 20
50-40922	Rigid	Single-Use Flat Brush	Light Blue	220 mm	60 mm x 20 mm x 15 mm	Bag of 20


Reusable models are also available in bags of 50. Add the suffix "- 50".

BATRIK

REUSABLE ORTHOPEDIC RASPS CLEANING BRUSH

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



PRODUCT DESCRIPTION

ITEM #	TOTAL LENGTH	HANDLE LENGTH	BRUSH SIZE	PACKAGING
50-RASP	9" / 23 cm	5.5" / 14 cm	0.78" x 2.7" x 3.5" (2 cm x 7 cm x 9 cm)	3 / Bag


BATRIK

REUSABLE PRECISION BRUSH

INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES
THE NEED FOR THE STAINLESS
STEEL TOOTHBRUSH-STYLE BRUSH,
WHICH CAN DAMAGE VALUABLE
INSTRUMENTS!

PRODUCT SELECTION

BLUE DELRIN PLASTIC

ITEM #	DESCRIPTION	TOTAL LENGTH
50-40800	Reusable Precision External Cleaning Brush	6.29" / 160 mm

PACKAGING: Bag of 3



BATRIK


REUSABLE INSTRUMENT CLEANING BRUSHES

Designed specifically for hospital applications, our Reusable Instrument Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

REUSABLE JAR/BOTTLE BRUSH

STERILIZER DRAIN BRUSH



PRODUCT SELECTION

BLUE DELRIN PLASTIC

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH
50-10295	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	7 mm ± 2 mm	13.19" / 355 mm
50-10224	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm

PACKAGING: Bag of 3

BATRIK ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.



STRAIGHT BRUSH



ROUND BRUSH



PRODUCT SELECTION

GREEN DELRIN PLASTIC

ITEM #	DESCRIPTION	HANDLE LENGTH	STRAIGHT BRUSH	ROUND BRUSH HEAD
50-ACETAB	Acetabular Reamer Cleaning Brush	4.72" / 120 mm ø 0.47" / 12 mm	1.18" / 30 mm ø 0.3125" / 8 mm	2" / 50 mm ø 1.75" / 45 mm

PACKAGING: Bag of 3

NOTE: Specifications subject to small variances

BATRIK

ALL PURPOSE REUSABLE BRUSHES

Plastic Handle Brushes

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

In-Brush™ Plastic Handle

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.

**28-72550****28-69394****29-72550****29-69394****29-70050****50-40911**

PRODUCT SELECTION

ITEM #	DESCRIPTION	LENGTH TOTAL	BRUSH LENGTH	BRUSH HEIGHT
28-72550	Plastic Single End Nylon Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
29-72550	Plastic Single End Stainless Steel Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
29-70050	Plastic Single End Brass Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
28-69394	Plastic Double End Nylon Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
29-69394	Plastic Double End Stainless Steel Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
50-40911	Hand Scrub Nylon Brush Reusable autoclavable	4.09" / 104 mm	3.74" / 95 mm	1.1" / 28 mm

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.

BATRIK

AUTOCLAVE CLEANING BRUSH KITS



ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

PRODUCT SELECTION

REGULAR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVEREG	Autoclave REGULAR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8" / 8.89 cm x 20.32 cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm

SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVESUP	Autoclave SUPERIOR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8" / 11.5 cm x 20.32cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm
		SF-100	2 SAVEFACE™ Face Shields	9" / 23 cm

REPLACEMENT COMPONENTS



HAND BRUSH
50-35900

Packaging: 1



BRUSH HEAD
50-55800

Packaging: 1



EXTENSION POLE
50-EXTPOLE

Packaging: 1



SAVEFACE™
FACE SHIELDS SF-100

Packaging: 100 / Box

100

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

SURGMED

GROUP

Proudly Canadian, **The Surgmed™ Group of Companies** began operations in 1992 as a single distribution company known as **Advanced Surgi-Pharm Inc.** The company quickly became known as an industry leader and specialty distributor of value-added surgical products in the Canadian hospital market. Five years later, **Surgmed expanded upon its success with the launch of Batrik Medical Manufacturing Inc.** in 1997 and the subsequent acquisition of **Imperial Surgical Ltd.** one year later in 1998.

Together, our three unique companies serve the Canadian and global healthcare markets with high quality and competitively priced products, all backed by world-class customer care.



Founded in 1992, Advanced Surgi-Pharm Inc. is the distribution arm of Surgmed™ in Canada, representing leading US, international and domestic manufacturers for a wide range of specialty products with unique benefits for applications in the Operating Room, Women's Health, Labor & Delivery, ICU, and Sterile Processing areas. Equipped with a direct sales team, ASP provides effective nationwide coverage as the Canadian distributor of Batrik, Vygon, Rocket Medical, Cooper Surgical, SunMed, and other industry leading products including a wide-ranging portfolio of top-rated PPE products.

Established in 1997 and sold in 80 countries globally, Batrik Medical Manufacturing Inc. specializes in the manufacturing of endoscopy accessories, enzymatic solutions, infection control, central sterilization, and surgical devices. Batrik's in-house controlled environment production facility adheres to strictly enforced quality and regulatory control systems. They are registered as a device manufacturing establishment with the United States FDA and are in full compliance with CE and Health Canada Regulations. All Batrik products are rigorously tested, liability insured, and meet/exceed high international quality standards.

Established in 1935, Imperial Surgical Ltd. is one of Canada's oldest hospital capital equipment manufacturers, well-known for their stainless steel healthcare equipment for Sterile Processing, Operating Room, Examination and Treatment, Patient Handling, and Mortuary/ Pathology areas. Their growing portfolio of products also has applications in Infection Control & Prevention, an important area of focus for their future.



SURGMED GROUP

2165 46th Ave, Lachine, QC, Canada, H8T 2P1

Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax.: +1.514.631.9083

Email: info@surgmed.com | Website: www.batrik.com